

Product datasheet for **RC225900L3V**

GTF3C5 (NM_001122823) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GTF3C5 (NM_001122823) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GTF3C5
Synonyms:	TFIIC2-63; TFIIC63; TFIICepsilon
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001122823
ORF Size:	1578 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225900).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001122823.1
RefSeq ORF:	1581 bp
Locus ID:	9328
UniProt ID:	Q9Y5Q8
Cytogenetics:	9q34.13
Protein Families:	Transcription Factors
MW:	60.2 kDa



[View online »](#)

Gene Summary:

Involved in RNA polymerase III-mediated transcription. Integral, tightly associated component of the DNA-binding TFIIC2 subcomplex that directly binds tRNA and virus-associated RNA promoters.[UniProtKB/Swiss-Prot Function]