

## Product datasheet for **RC225906L4V**

### **ANKRD53 (NM\_001115116) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

|                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Symbol:</b>                   | ANKRD53                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Mammalian Cell Selection:</b> | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Vector:</b>                   | pLenti-C-mGFP-P2A-Puro (PS100093)                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Tag:</b>                      | mGFP                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>ACCN:</b>                     | NM_001115116                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ORF Size:</b>                 | 1590 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>ORF Nucleotide Sequence:</b>  | The ORF insert of this clone is exactly the same as(RC225906).                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>OTI Disclaimer:</b>           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| <b>OTI Annotation:</b>           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| <b>RefSeq:</b>                   | <a href="#">NM_001115116.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>               | 1593 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>                 | 79998                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>UniProt ID:</b>               | <a href="#">Q8N9V6</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Cytogenetics:</b>             | 2p13.3                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>MW:</b>                       | 59.4 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Gene Summary:</b>             | Required for normal progression through mitosis. Involved in chromosome alignment and cytokinesis via regulation of microtubules polymerization.[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                                                                                    |



