

Product datasheet for **RC225907**

USP17 (NM_001105662) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	USP17
Synonyms:	MGC119330; MGC119331; MGC119333; RS447; USP17A; USP17H; USP17I; USP17J; USP17K; USP17L; USP17M
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC225907 representing NM_001105662
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGGACGACTCACTCTACTTGGGAGGTGAGTGGCAGTTCAACCACTTTCCAAAACCTCACATCTTCTC
 GGCCCCGATGCAGCTTTTGTGAAATCCAGCGGACTTCTCTCCCTGAGAAGTCACCACTCTCATGTGAGAC
 CCGTGTGCGACTCTGTGACGATTTGGCTCTGTGGCAAGACAGCTTGTCTCCAGGGAGAAGCTTCTCTGT
 AGTAGCAGGAGACCTGCTGCGGTGGGGCTGGGCTCCAGAATATGGGAAATACCTGCTACGTGAACGCTT
 CCTTGGAGTGGGTGACATACACACCGCCCTTGCCAACTACATGCTGTCCCGGAGCACTCTCAAACGTG
 TCATCGTCAACAAGGGCTGCATGCTCTGTACTATGCAAGCTCACATCACACGGGCCCTCCACAATCTGGG
 CACGTATCCAGCCCTCACAGGCATTGGCTGCTGGCTTCCATAGAGGCAAGCAGGAAGATGCCCATGAAT
 TTCTCATGTTCACTGTGGATGCCATGAAAAGGCATGCCTTCCGGGCACAAGCAGGTAGATCATCACTC
 TAAGGACACCACCTCATCCACCAATATTTGGAGGCTACTGGAGATCTCAAATCAAGTGTCTCCACTGC
 CACGGCATTTCAGACACTTTTGACCTTACCTGGACATCGCCCTGGATATCCAGGCAGCTCAGAGTGTCC
 AGCAAGCTTTGGAACAGTTGGTGAAGCCGAAGAACTCAATGGAGAGAATGCCTATCATTGTGGTGTGTTG
 TCTCCAGAGGGCGCCGCTCCAAGACGTTAACTTTACACACCTCTGCCAAGGTCTCATCTTGTATTG
 AAGAGATTCTCCGATGTACAGGCAACAAGATTGCCAAGAATGTGCAATATCCTGAGTGCCTTGACATGC
 AGCCATACATGTCTCAGCAGAACACAGGACCTTGTCTATGTCCTCTATGCTGTGCTGGTCCACGCTGG
 GTGGAGTTGTCAACGACATTACTTCTTATGTCAAAGCTCAAGAAGGCCAGTGGTATAAAATGGAT
 GATGCCGAGGTACCGCTCTAGCATCACTTCTGTCTGAGTCAACAGGCCTACGTCTCTTTTACATCC
 AGAAGAGTGAATGGAAAGACACAGTGAGAGTGTGCAAGAGCAGGGAACCAAGAGCCCTTGGCGCAGA
 AGACACAGACAGGCGAGCAACGCAAGGAGAGCTCAAGAGAGACCACCCCTGCCTCCAGGCCCCCGAGTTG
 GACGAGCACTTGGTGGAAAGAGCCACTCAGGAAAGCACCTTAGACCACTGGAATTCCTTCAAGAGCAAA
 ACAAAACGAAGCCTGAGTTCAACGTGAGAAAGTCAAGGTACCTGCCTCCGACGTACTTGTGATTCA
 TCAATCAAAATACAAGTGTGGATGAAGAACCATCATCCTGAACAGCAAGCTCCCTGCTAAACCTCTCT
 TCGACGACCCCGACACATCAGGAGTCCATGAACACTGGCACACTCGCTTCCCTGCGAGGGAGGCCAGGA
 GATCCAAGGGAAGAACAACACAGCAAGAGGGCTCTGCTTGTGTGCCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225907 representing NM_001105662
 Red=Cloning site Green=Tags(s)


MEDDSLVLGGWEQFNHFKLTSSRPDAFAEIQRSLPEKSPLSCETRVDLCDDLPVARQLAPREKLPL
 SSRPRAAVGAGLQNMGNTCYVNASLEWLTYPPLANYMLSREHSQTCHRHKGCMCTMQAHITRALHNPG
 HVIQPSQALAAGFHRGQEDAHEFLMFTVDAMKKAACLPGHKQVDHHSKDTTLIHQIFGGYWRSQIKLHC
 HGISDTFDPYLDIALDIQAASVQQAELVKPEELNGENAYHCGVCLQRAPASKTLTLHTSAKVLILVL
 KRFSVDVTGNKIAKNVQYPECLDMQPYMSQQNTGPLVYVLVAVLVHAGWSCHNGHYFSYVKAQEGQWYKMD
 DAEVTASSITSVLQQAYVLFYIQKSEWERHSESVSRGREPRALGAEDTDRRATQGELKRDHPCLQAPEL
 DEHLVERATQESTLDHWKFLQEQNKTKPEFNVRKVEGTLPPDVLVIHQSKYKCGMKNHHPQQSSLLNLS
 STTPTHQESMNTGTLASLRGRARRSKGNKHSKRALLVCQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8035_a08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001105662

ORF Size: 1590 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001105662.1](#), [NP_001099132.1](#)


RefSeq ORF: 1592 bp

Locus ID: 391627
Cytogenetics: 4p16.1
MW: 59.6 kDa

Product images:



Circular map for RC225907



Western blot validation of overexpression lysate (Cat# [LY426293]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225907 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).