

Product datasheet for **RC225912**

ZNF837 (NM_138466) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ZNF837
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC225912 representing NM_138466
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGGCTCCAGCCCAGAAGGCTGGGCAGGGAGGACTCCCCAAGGCCGATGCCAGGGTGCCTCCGGGG
 CCCGGGAGAAGAGGCCCGAGGAGCCGAGGCCCTCGAAGAGGACCGAGCTGGGAGCCGCCCACTCAAAA
 GGGGGACCTGCGTGGAGCGCGGGCGGTAGGACGACTCCCCAGGCGGGGGCTCCCGGGGTGCAGCCTC
 GGGGTGAGCCCCGGCCCGGGGACCCGGCACAGCGCCGGGACCAGACCCCTCGTGCGGAGCCGTGCGGCC
 CCACCTCTTCTCAGAACCCTGAGCTGGTTATCCCCGAGGGGCTGCAAGCCGGGAGGGCCCTGTAGGAG
 CCCAGGCGGTGGCGGGGACTGCAGCAGGAACTCCTGCCTGGCGTGGCATCGCGGGGCGCCCGCTGGGGAG
 ACGCCACCCGTGTGTGACCCCTGTCCGGAGCGGATCCAGAACCACCCCGGACTCAACTGTGTGAGGTCC
 ACACGACTGTTGGCGTGCCAAACAGGGACTGGCGCTCCGACCTGCCGAGGACCCAAAGCCGACCTC
 CCGCGGAGGAACCCCTTGGTGGAGCAGCCCGGGCTTGGCGTGGCGGAGGCCTTTGCGTGGAGGGCC
 CTGCGGATCCCCAGGAGAGGCTGCAGGCGACGAGGAGCCCGTCCGTGTGCCCGGTGCGGGAAGCGCT
 TCCGCCCAACAGCAGCAGAGGCGGGCAAGAGTCCCCAGTGTGTCTGAGTGGCGCAAACCTCGCG
 ACCTCGCCCTATTGTCCCGACCCCGGGCCAGCGACTGTACGCTTGCAGCAGTGGCGCAAGGCCTTC
 ACGCGACCTCCAGCCTGCTGCAGCACCAGCGCATCCACACGGGCGAGCGGCCCTACGAGTGCGCCGAGT
 GCGGCAAGGCCTTGTGCGCTGCTCCGGCCTGTACCGCCACCAAGACGCACTCGGCCGAGCGCCACCG
 GCGTGGCCCCGTCTGCGCGGCGCGCCTTCCGGCTGGGGTGGCCGCCCTGCGGGGACTACAGCAGCGG
 AGTCCCCGACGGGGTGGGAGCCGGGGAGAAGCCGTACGAGTGCGCCGACTGCGCCAAGGCCTTCGGGC
 TGTTCTCGACCTCGTGGAGCACCAGCGCGTGCACACCGCGGAGAAGCCCTACGCGTCCCCGAGTGGCG
 CAAGGCCTTCAACCAGCGCTCGAACCCTGAGCCGGCACCAGCGCACTCACAGCAGCGCCAAGCCCTACGCG
 TGCCCACTGTGCGAAAAGGCCTTCAAGGGCCGCTCGGGCCTGGTGAACACCAGCGCGCGCACACCGGCG
 AGCGGCCCTATGGTGTCTCCGAGTGCAGAAAGACCTTCCCGGGTGTCTCCGAGTGCGCCAGCAGAGCG
 CCTGCACTCGGGCGAGAAGCCCTACATCTGCCGCGACTGCGGCAAGGCCTTGTGCGCAACTGCAGCCTG
 GTGCGCCACCTGCGCACGCACACGGGCGAGCGGCCCTACGCGTGGGAGATTGCGGCCGCGCCTTCAGCC
 AACGCTCCAACCTCAACGAGCACCGAAGCGGCACGGGGGCGCGCCGCGCCT

ACGCGACCGACGCTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225912 representing NM_138466
 Red=Cloning site Green=Tags(s)


MEAPAQKAGQGGLPKADAQGASGAREKRPEEPRPLEEDRAGSRPTQKGDLRGAAGGRTPPGGSGRCSL
 GVSPGPGTRHSAGTRPLVREPCGPTSSQNPELVIPEGLQAREGPCRSPARGGDCSRNSCLAWHRGAPAGE
 TPPVCDPCPERIQNHPRQLCEVHTDCWPCQPGTGAPTCPRTPKPTRSRGRNPLVEQPRACACGEAFAWRA
 LRIPQERLQATEEPRPCARCGKRFNPNQQQAGKSPPVCECGQTSRPRPIVPDPPAQRLYACDECGKAF
 TRTSSLQHQRHITGERPYEACGKAFVRCGLYRHKTHSAERHRRGPVLARRAFRLGCPPCGDYSER
 SPRRGSGAGEKPYECADCAKAFGLFSLVEHRRVHTGEKPYACPECGKAFNQRSNLSRHQRTHSSAKPYA
 CPLCEKAFKGRSLVQHQAHTGERPYGCSECGKTFRGCELRQHERLHSGEKPYICRDCGKAFVRNCSL
 VRHLRTHTGERPYACGDCGRAFSQRSNLSNEHRKRHGGRAAP

SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8104_h09.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_138466

ORF Size: 1593 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_138466.2](#)


RefSeq Size: 1948 bp

RefSeq ORF:	1596 bp
Locus ID:	116412
UniProt ID:	Q96EG3
Cytogenetics:	19q13.43
MW:	58.1 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC225912



Western blot validation of overexpression lysate (Cat# [LY430049]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC225912 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).