

Product datasheet for **RC225919L3V**

FAM196B (INSYN2B) (NM_001129891) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	INSYN2B
Synonyms:	C5orf57; FAM196B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001129891
ORF Size:	1605 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225919).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001129891.1
RefSeq ORF:	1608 bp
Locus ID:	100131897
UniProt ID:	A6NMK8
Cytogenetics:	5q35.1
MW:	59 kDa

