

## Product datasheet for **RC225954L3V**

### **TBCCD1 (NM\_001134415) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TBCCD1 (NM_001134415) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TBCCD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001134415
ORF Size:	1671 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225954).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001134415.1</a> , <a href="#">NP_001127887.1</a>
RefSeq Size:	3080 bp
RefSeq ORF:	1674 bp
Locus ID:	55171
UniProt ID:	<a href="#">Q9NVR7</a>
Cytogenetics:	3q27.3
MW:	63.5 kDa
Gene Summary:	Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.[UniProtKB/Swiss-Prot Function]



[View online »](#)