

Product datasheet for RC225967L3V

NEDD1 (NM_001135177) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NEDD1 (NM_001135177) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEDD1
Synonyms:	GCP-WD; TUBGCP7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135177
ORF Size:	1983 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225967).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001135177.1 , NP_001128649.1
RefSeq Size:	3738 bp
RefSeq ORF:	1716 bp
Locus ID:	121441
UniProt ID:	Q8NHV4
Cytogenetics:	12q23.1
MW:	72 kDa


[View online »](#)

Gene Summary: Required for mitosis progression. Promotes the nucleation of microtubules from the spindle.
[UniProtKB/Swiss-Prot Function]