

Product datasheet for **RC225993**

FIPIL1 (NM_001134937) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FIPIL1
Synonyms:	FIP1; hFip1; Rhe
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RC225993 representing NM_001134937
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTCGGCCGCGAGGTCGAGCGCCTAGTGTGCGAGCTGAGCGCGGGACCGAGGGGATGAGGAGGAAG
 AGTGGCTCTATGGCGATGAAAATGAAGTTGAAAGGCCAGAAGAAGAAATGCCAGTGCATACCTCCATC
 TGGAAATTGAAGATGAACTGCTGAAAATGGTGTACAAAACCGAAAGTGACTGAGACCGAAGATGATAGT
 GATAGTGACAGCGATGATGATGAAGATGATGTTTCATGTCACATATAGGAGACATTAACCGGGAGCACCAC
 AGTATGGGAGTTATGGTACAGCACCTGTAATCTTAACATCAAGACAGGGGAAGAGTTTATGGAAGTAC
 AGGGACAAAAGTCAAAGGAGTAGACCTTGATGCACCTGGAAGCATTAAATGGAGTTCCACTCTTAGAGGTA
 GATTTGGATTCTTTTGAAGATAAACCATGGCGTAAACCTGGTGTGATCTTTCTGATTATTTTAATTATG
 GGTTTAAATGAAGATACCTGGAAGCTTACTGTGAAAAACAAAAGAGGATACGAATGGGACTTGAAGTTAT
 ACCAGTAACCTCTACTACAAAATAAATTACGGCCGAAGACTGTACTATGGAAGTTACACCAGGTGCAGAG
 ATCCAAGATGGCAGATTCAATCTTTTAAAGGTACAGCAGGGAAGAACTGGAAGTACAGAGAAAGAACTG
 CCCTTCCATCTACAAAAGCTGAGTTTACTTCTCCTCCTCTTTGTTCAAGACTGGGCTTCCACCGAGCAG
 AAACAGCACTTCTTCTCAGTCTCAGACAAGTACTGCCTCCAGAAAAGCCAATTAACGCGTTGGGAAGTGG
 CAGGATCGATATGGGAGGGCCGAATCACCTGATCTAAGGAGATTACCTGGGCAATTGATGTTATCGGTC
 AGACTATAACTATCAGCCGAGTAGAAGGCAGGCGACGGGCAATGAGAACAGCAACATACAGGTCTCTTTC
 TGAAAGATCTGCTACTGAAGTAGACAACAATTTAGCAAAACCACTCCGTTTTTCCCTCCAGGAGCTCCT
 CCCACTCACCTCCACCTCCTCCATTTCTCCACCTCCTCGACTGTCAGCACTGCTCCACCTCTGATTC
 CACCACCGGATTTCCAATACTGTACCACCTCCAGGTTTTCTCCTCCACCGGCTCCACCTCCATC
 TCTTATACCAACAATAGAAAGTGGACATTCTCTGGTTATGATAGTCGTTCTGCACGTGCATTTCCATAT
 GGCAATGTTGCCTTTCCCATCTTCTCTGTTCTGCTCCTTCGTGGCCTAGTCTTGTGGACACCAGCAAGC
 AGTGGGACTATTATGCCAGAAGAGAGAAAGACCGAGATAGAGAGAGAGACAGAGACAGAGAGCGAGACCG
 TGATCGGGACAGAGAAAGAGAACGCCACAGAGAGAGAGAGGGAGCGTGATCACAGTCCTACACCAAGT
 GTTTTCAACAGCGATGAAGAACGATACAGATACAGGGAATATGCAGAAAGAGGTTATGAGCGTCACAGAG
 CAAGTCGAGAAAAAGAAGAACGACATAGAGAAAGACGACACAGGGAGAAAGAGGAAACAGACATAAGTC
 TTCTCGAAGTAATAGTAGACGTGCGCATGAAAGTGAAGAAGGAGATAGTCACAGGAGACACAACACAAA
 AAATCTAAAAGAAGCAAGAAGGAAAAGAAGCGGCAGTGAGCCTGCCCTGAACAGGAGAGCACCGAAG
 CTACACCTGCAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225993 representing NM_001134937
 Red=Cloning site Green=Tags(s)

MSAGEVERLVSEL SGGTGGDEEEWLYGDENEVERPEEENASANPPSGIEDETAENGVPKPKVTETEDDS
 DSDSDDDEDDVHVTIGDIKTGAPQYGSYGTAPVNLNIKTGGRVYGTGTGKVKGVLDLAPGSINGVPLLEV
 DLDSFEDKPWRKPGADLSDYFNYGFNEDTWKAYCEKQKRIRMGLEVIPVTSTTNKITAEDCTMEVTPGAE
 IQDGRFNLKVKVQGRGTGNSEKETALPSTKAEFTSPPSLFTGLPPSRNSTSSQSQTSTASRKANSSVGKW
 QDRYGRAESPDRLRLPGAIDVIGQITITISRVGRRRANENSNIQVLSERSATEVDNNFSKPPPFPPPGAP
 PTHLPPPPFLPPPTVSTAPPLIPPGIPITVPPPGFPPPGAPPPSLIPTIESGHSSGYDSRSARAFPY
 GNVAFPHLPGSAPSWPSLVDTSKQWDYYARREKDRDRERDRDRERDRERERTRERERERDRHSPTPS
 VFNSDEERYRYREYAERGYERHRASREKEERHRERHREKEETRHKSSRSNSRRRHESEEGDSHRRHKHK
 KSKRSKEGKEAGSEPAPEQESTETPAE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8034_c02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001134937

ORF Size: 1764 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001134937.2](#)

RefSeq ORF: 1767 bp

Locus ID: 81608


UniProt ID: [Q6UN15](#)

Cytogenetics: 4q12

MW: 65.5 kDa

Gene Summary: This gene encodes a subunit of the CPSF (cleavage and polyadenylation specificity factor) complex that polyadenylates the 3' end of mRNA precursors. This gene, the homolog of yeast Fip1 (factor interacting with PAP), binds to U-rich sequences of pre-mRNA and stimulates poly(A) polymerase activity. Its N-terminus contains a PAP-binding site and its C-terminus an RNA-binding domain. An interstitial chromosomal deletion on 4q12 creates an in-frame fusion of human genes FIP1L1 and PDGFRA (platelet-derived growth factor receptor, alpha). The FIP1L1-PDGFRA fusion gene encodes a constitutively activated tyrosine kinase that joins the first 233 amino acids of FIP1L1 to the last 523 amino acids of PDGFRA. This gene fusion and chromosomal deletion is the cause of some forms of idiopathic hypereosinophilic syndrome (HES). This syndrome, recently reclassified as chronic eosinophilic leukemia (CEL), is responsive to treatment with tyrosine kinase inhibitors. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008]

Product images:



Circular map for RC225993