

Product datasheet for **RC226009**

PEX5 (NM_001131024) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX5 (NM_001131024) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PEX5
Synonyms:	PBD2A; PBD2B; PTS1-BP; PTS1R; PXR1; RCDP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC226009 representing NM_001131024
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCAATGCGGGAGCTGGTGGAGGCCAATGCGGGGTGCCAACCCGCTCATGAAGCTCGCCGGGCACT
 TCACCCAGGACAAGGCCCTTCGGCAGGAGGATTGAGGCCTGGCCCTGGCCCCGGAGCCCCGGCCTC
 TGAGGCAGCCTCCAAGCCTTTGGGAGTAGCTTCTGAAGATGAGTTGGTGGCTGAATTCCTGCAGGACCAG
 AATGCACCCCTTGTGTCCCGTGCCTCAGACCTCAAGATGGATGACCTCCTGGCTGAGATGCAGCAGA
 TTGAGCAGTCAAACCTCCGCCAGGCTCCCAGAGAGCCCTGGTGTGGCAGACTTGGCCTTGTCTGAGAA
 CTGGGCCAGGAGTTTCTGCAGCTGGAGATGCTGTGGATGTAACCTCAGGATTATAATGAGACTGACTGG
 TCCCAAGAATTCATCTCTGAAGTTACAGACCCCTTGTCTGTGTCCCTGCCCGTGGGCTGAGGAATATT
 TGGAGCAATCAGAGGAGAAGCTGTGGCTGGGAGAACCTGAGGGAACAGCCACCGATCGCTGGTATGATGA
 ATATCATCTGAGGAGGATCTGCAGCACACGGCCAGTGACTTTGTGGCCAAAGTGGATGACCCCAAATTG
 GCTAATTCGAGGGTACATCAGATGCCTGGGTTGACCAGTTCACAAGACCAGTAAACACATCTGCCCTTG
 ATATGGAGTTTGAACGAGCCAAGTCAGCTATAGAGTCTGATGTCGATTTCTGGGACAAGTTGCAGGCAGA
 GTTGGAGGAGATGGCAAACGGGATGCTGAGGCCACCCCTGGCTTTCTGACTATGATGACCTTACGTCA
 GCTACCTATGATAAGGGGTACCAGTTTGGAGGAGAAACCCCTTGCCTGATCACCTCAGCCTTTTGAAG
 AAGGGCTGCGGCGCCTTCAGGAGGGGACTGCCAAATGCTGTGTGCTTTTGGAGCAGCTGTGCAGCA
 GGATCCTAAGCACATGGAAGCTTGGCAGTATCTGGGTACCACCCAGGCAGAGAATGAACAAGAATTA
 GCCATCAGTGCATTGCGGAGGTGTCTGGAGCTAAAGCCAGATAACCAGACAGCACTGATGGCGCTGGCTG
 TGAGCTTACCAACGAGTCCCTGCAGCGACAGGCTGTGAAACCCACGAGACTGGCTGCGGTACACACC
 AGCCTATGCCCATCTGGTGACACCTGCTGAAGAAGGGGCTGGTGGGCGAGGACTGGGCCCCAGCAAGCT
 ATCCTGGGATCTCTTGTCTGACTCCCTGTTTCTGAAGTAAAAGAGCTCTTCTGGCAGCTGTGCGGC
 TGGACCTACCTCCATTGACCCTGATGTGCAGTGTGGCTTGGGAGTCTTTTCAACCTGAGTGGGGAGTA
 TGACAAGGCCGTGGACTGCTTACAGCTGCCCTCAGCGTTCGTCCTCAATGACTATTTGCTGTGGAATAAG
 CTAGGCGCCACCCTGGCCAATGGAACCAGAGTGAAGAAGCAGTACGTCGTACCGCCGGGCCCTCGAGC
 TCCAGCCTGGCTATATCCGGTCCCGCTATAACCTGGGCATCAGCTGCATCAACCTCGGGGCTCACCGGA
 GGCTGTGGAGCACTTTCTGGAGGCCCTGAACATGCAGAGGAAAAGCCGGGGCCCCGGGGTGAAGGAGGT
 GCCATGTCGGAGAACATCTGGAGCACCTGCGTTTGGCATTGTCTATGTTAGGCCAGAGCGATGCCTATG
 GGCAGCCGACGCGGGGATCTGTCCACCTCCTAACTATGTTTGGCCTGCCCCAG

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC226009 representing NM_001131024
 Red=Cloning site Green=Tags(s)

MAMRELVEAECGGANPLMKLAGHFTQDKALRQEGLRPGWPPGAPASEAASKPLGVASEDELVAEFLQDQ
 NAPLVSRAPQTFKMDLLAEMQQIEQSNFRQAPQRAPGVADLALSENWAQEF LAAGDAVDVTDQYNETDW
 SQEFLISEVTDPLSVSPARWAEYLEQSEKWLWGEPEGTATDRWYDEYHPEEDLQHTASDFVAKVDDPKL
 ANSEGTSDAWVDQFTRPVNTSALDMEFERAKSAIESDVDFWDLQAELEEMAKRDAEAHPWLSYDLDLTS
 ATYDKGYQFEENPLRDHPQPFEEGLRRLQEGDLPNAVLLFEAAVQDDPKHMEAWQYLGTTQAENEQELL
 AISALRRCLELKPNDQALMALAVSFTNESLQRQACETLRDWLRYTPAYAHLVTPAEEGAGGAGLGP SKR
 ILGSLSDSLFLEVKELFLAAVRLDPTSIDPDVQCGLGVLFNLSGEYDKAVDCF TAALSVRPNDYLLWNK
 LGATLANGNQSEEAVAA YRRALELQPGYIRSRYNLGISCNLGAHREAVEHFLEALNMQRKSRGPRGEGG
 AMSENIWSTLRLALSMLGQSDAYGAADARDLSTLLTMFGLPQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8069_c06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001131024

ORF Size: 1806 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001131024.2](#)

RefSeq ORF: 1809 bp

Locus ID: 5830

UniProt ID: [P50542](#)


Cytogenetics: 12p13.31

Protein Families: Druggable Genome


MW: 66.6 kDa

Gene Summary: The product of this gene binds to the C-terminal PTS1-type tripeptide peroxisomal targeting signal (SKL-type) and plays an essential role in peroxisomal protein import. Peroxins (PEXs) are proteins that are essential for the assembly of functional peroxisomes. The peroxisome biogenesis disorders (PBDs) are a group of genetically heterogeneous autosomal recessive, lethal diseases characterized by multiple defects in peroxisome function. The peroxisomal biogenesis disorders are a heterogeneous group with at least 14 complementation groups and with more than 1 phenotype being observed in cases falling into particular complementation groups. Although the clinical features of PBD patients vary, cells from all PBD patients exhibit a defect in the import of one or more classes of peroxisomal matrix proteins into the organelle. Defects in this gene are a cause of neonatal adrenoleukodystrophy (NALD), a cause of Zellweger syndrome (ZWS) as well as may be a cause of infantile Refsum disease (IRD). Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]

Product images:



Circular map for RC226009



Western blot validation of overexpression lysate (Cat# [LY427362]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226009 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).