

Product datasheet for **RC226021**

Phostensin (PPP1R18) (NM_001134870) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Phostensin (PPP1R18) (NM_001134870) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Phostensin
Synonyms:	HKMT1098; KIAA1949
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC226021 representing NM_001134870
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGCCACCATCCCAGACTGGAAGCTACAGCTGCTAGCCCGCGCCGCGCAGGAGGAGCGTCCGTTTCGAG
 GCCGAGAGAAAGCAGAACGGGAGCGCCTGTCCAGATGCCAGCCTGGAACGAGGGCTCCTGGAGCGCCG
 CCGGGCCAAGCTTGGGCTGTCCCTGGGAGCCTAGCCCTGTGCTAGGGACTGTAGAGGCTGGACCTCCA
 GACCCGGATGAGTCTGCGGTCTTCTGGAGGCCATCGGGCAGTGCACCAGAACCATTATCCGGCAGG
 AGCGGCAGCAGCAGCAGCAACAACCGAGTGAAGAGCTGTAGCAGAGAGAAAGCCTGGGCTCT
 GGAGGCCCGGGAGCGGAGACCCAGCCCTGGGAGATGCGGGATCAGAGCCCCAAGGAAGAGAGTCAAGA
 GAAGAGAGACTAAGTCCGAGGGAGACCAGAGAGAGGAGGCTGGGATAGGGGAGCCCAAGATTGAGCC
 TGAGGCCCTCTGGAGGCTCGGACTGGAGGCAAGCCAGGAGAGGTGGGAGACAGGAGCTCCCGACTGTC
 AGAGGCATGGAATGGAGGCTGAGTCTGGAGAACTCCAGAGCGGAGTCTGAGACTAGCAGAGTCTCGA
 GAGCAAGCCCCAGGAGAAAAGAGGTGGAAGTAGACTGAGCCAGGGGAATCTGCCTACCAGAAGTTGG
 GCCTGACAGAGGCCCATAAATGGAGACCTGACTCCAGAGAGTCTCAGGAACAGAGTTTGGTACAACCTGGA
 GGCAACAGAGTGGAGGCTGAGGTGAGGAGAGAAAGACAAGACTACTCGGAAGAATGTGGGAGAAAAGAA
 GAGTGGCCAGTTCCAGGGGTAGCTCCAAAAGAGACTGCAGAGCTGTCCGAGACCTTGACAAGGGAGGCC
 AAGGCAACAGTTCCGAGGAGTGGAGGCAGCAGAGCAGAGGCCTGTGGAAGATGGCGAGAGGGGCATGAA
 GCCAACAGAAGGTGGAATGGACCTGAATCCGGGAAGGCTCGAGAATGGACCCAGGGACATAGAG
 GCTCAAACCTCAGAAACCAGAACCTCCAGAGTCAAGAGAGAGCTTCTGGAATCTCCCGGTGTGGAGGCTG
 GAGAAGGGGAGGCTGAGAAGGAGGAGCGGGGCTCAGGCAGGCCTCTGAGAGCCCTGCAGAACTGCTG
 CTCTGTGCCCTCCCCCTCCACACAGAGACGCTGGGACTGGAGGCTGAGACAGCAGGAAGAGGAAGCA
 GTGGAGTCCAGCCCCACACAGCCCTCTGTCTCCCCACCCCCAGCCCCAACTGCCCCCAACCTC
 CTGGGGATCCCCTCATGAGCCGCTGTTCTATGGGGTGAAGGCAGGGCCAGGGTGGGGGCCCCCGCCG
 CAGTGGACACACCTTACCGTCAACCCCGCGGTCTGTGCCCTGCGACCCAGCCACCCCAACCTCT
 CCAGCCACAGTTGATGCTGCAGTCCCGGGGCTGGGAAGAAGCGGTACCAACTGCCGAGGAGATCTTGG
 TTCTGGGGGCTACCTCCGTCTCAGCCGAGCTGCCTTGCCAAGGGTCCCCGAAAGACACCACAAACA
 GCTTAAGATCTCCTTCAAGCAGACAGCCCTGGAGACCAGTACCAATACCCCTCCGAGAGTTCGGTACTG
 GAGGAGCTGGGCCGAGCCTGAGGTCCCAAGTCCCCCAACCTCCAGAGCCCAACCCGACGACGAAG
 AGGATGAGGAAGAGCTGCTGCTGCTGCAGCCAGAGCTCCAGGGCGGGCTGCGACCAAGGCCTGATTGT
 GGATGAGTCTGCCGGCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC226021 representing NM_001134870
 Red=Cloning site Green=Tags(s)


MATIPDWKLQLLARRRQEEASVRGREKAERERLSQMPAWKRGLLERRRAKLGLSPGEPSPVLGTVEAGPP
 DPDESALLLEAIGPVHQNRFIQERQQQQQQRSEELLAERKPGPLEARERRPSGEMRDQSPKGRESR
 EERLSPRETRERRLGIGGAQELSLRPLEARDWRQSPGEVGRSSRLSEAWKWRLSPGETPERSLRLAESR
 EQSPRRKEVESRLSPGESAYQKLGLTEAHKWRPDSRESQEQSLVQLEATEWRLRSGEERQDYSEECGRKE
 EWPVPGVAPKETAESETLTREAQGNSSAGVEAAEQRPVEDGERGMKPTGKWTLSNGKAREWTPRDIE
 AQQTQKPEPPESAELLESPPGVEAGEGEAEKEEAGAQRPLRALQNCSSVPSPLPPEDAGTGGLRQEEEA
 VELQPPPPAPLSPPPAPTAPQPPGDPLMSRLFYGVKAGPGVGAPRRSGHTFTVNPRRSVPPATPATPTS
 PATVDAAVPGAGKKRYPTAEIILVGGYLRSLRSLAKGSPERHHKQLKISFSETALETYYQYPSESSVL
 EELGPEPEVPSAPNPAAQPDDEEDEEELLLQPELQGLRTKALIVDESCRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8028_e03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001134870

ORF Size: 1839 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001134870.2](#)

RefSeq Size: 2969 bp

RefSeq ORF: 1842 bp
 Locus ID: 170954
 UniProt ID: [Q6NYC8](#)
 Cytogenetics: 6p21.33
 MW: 67.9 kDa

Gene Summary: Protein phosphatase-1 (PP1; see MIM 176875) interacts with regulatory subunits that target the enzyme to different cellular locations and change its activity toward specific substrates. Phostensin is a regulatory subunit that targets PP1 to F-actin (see MIM 102610) cytoskeleton (Kao et al., 2007 [PubMed 17374523]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RC226021