

Product datasheet for **RC226025L3V**

SLC6A12 (NM_001122847) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC6A12 (NM_001122847) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC6A12
Synonyms:	BGT-1; BGT1; GAT2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001122847
ORF Size:	1842 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226025).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001122847.1
RefSeq Size:	3165 bp
RefSeq ORF:	1845 bp
Locus ID:	6539
UniProt ID:	P48065
Cytogenetics:	12p13.33
Protein Families:	Druggable Genome, Transmembrane
MW:	69.4 kDa



[View online »](#)

Gene Summary:

Transports betaine and GABA. May have a role in regulation of GABAergic transmission in the brain through the reuptake of GABA into presynaptic terminals, as well as in osmotic regulation.[UniProtKB/Swiss-Prot Function]