

Product datasheet for **RC226033L3V**

LRRC70 (NM_181506) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	LRRC70
Synonyms:	SLRN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181506
ORF Size:	1866 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226033).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_181506.4
RefSeq ORF:	1869 bp
Locus ID:	100130733
UniProt ID:	Q7Z2Q7
Cytogenetics:	5q12.1
MW:	70.1 kDa



Gene Summary:

Renders cells highly sensitive to the activation by cytokines and lipopolysaccharide (LPS).
[UniProtKB/Swiss-Prot Function]