

Product datasheet for **RC226040**

LRRIQ3 (NM_001105659) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	LRRIQ3
Synonyms:	LRRC44
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC226040 representing NM_001105659
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTTTCATGGAACAGTCACAGAAGAGCTAACCAGTCATGAAGAATGGAGTCACTATAATGAAAACATAA
 GAGAAGGTCAAAAAGATTTTGTGTTTGTGAAGTTCAATGGCCTTCATTTAAAGTCTATGGAAAATTTGCA
 GTCTTGATCTCTCTTAGAGTATGCATCTTCTCAACAATTTTATTACAGATATTATCCACTTCAAAGT
 TGTATAAAATTAATCAAACTTGATCTCCATGGAATCAGATAAAGAGTCTACCAATACCAAAATTTTGA
 ATGGATTGAAGAACCTAACTACTCTATCTTCATGACAATGGGTTTGCAAAGTTAAAGAATATATGTGT
 ATTATCTGCCTGTCCAACCTCATTGCCCTCACTATGTTGATTGTCCAGTAAGCCTTAAAAAGGATAT
 AGACATGTTCTTGTAAACAGTATATGGCCTCTCAAGCGCTGGATCATCATGTGATTCTGATGAAGAAA
 TAATTCAGAACTGGCATCTTCTGAAAGATTCAAAGCATGTAACCATCGACTTTTCTTAATTTCTGCC
 AGCTTTGAGAAAGGGAACAACCTATGAAGAGGAAATTAATAATTTAAACATATTACTTCAAAAATTAAT
 GCAATTCTGGCTCATAATTACCAAGTTTGTATTGTTCAAAGATGGATACGTGGTTTCTTAGTTAGAAAAA
 ATTTGAGCCCTGTGTTTTCCACAAAAAAAACAGCAGGAAAAAATTATTAGAGGATATGAAGCAAAATG
 GATTTACATAACCAAGGGTATGAAGATAAGCTCCTTAAGGATCTCTTTTCAAACCTGAACTAATATA
 AAAGGAAAGCTTGATATTGGAACATAATATATATTATCCTGTTGATTAAAAATTCAGTGAACACA
 GAAAAACATGTGTATCTATTTTATGTGAATTAACCAAAAGATCTAGGCATGAAATCAAAACATCAAG
 ACATCTCATTCAAAAAGGTGAGGAGTCTGAAGATGAAATTGTGGATGAAAAATTGGATACCAAGTTTAGG
 ATATCAGTTTTCAAACCTACCATATACACTTCAGGTTCAATGAAGAATAATGCAGTATTGAGAGAGAAAA
 AACAACATTTTTCTGCATATCTCAGCCAATCTATACTACTCATCAAAGCCAATCATTAAAAAGA
 CATACGATTGGAGCGGAGTATGAAAGAGTTTTTGCACCACAAAGAGCTGGTATGAAACTCCGAACATTT
 AGTGATATTGACAAATATTACACAGAACAAAAGAAGCAGGAATACCATAAGAAAAAGTAAGAGTTGTAG
 CCATGGCACAAGTTGCTCGAGAAAGAGTTAGAGTAGCTGTTAATGAACATTTGAATCAGAAAAATATGC
 TACACAAAACTAATTGAAGAAAAAAGAGACAATTCAGAACAGTTTACGACAAGTTTGGCAAAACAGA
 TTCAACTACCTTGAAAAAGCTAGAGAGAGAAAAAGCTCTGTTTCTAAAAAGAAAAATCCCAAAAGGCTTCA
 AGCGTTTATTAGTTCAAACTTAAATAATGAGCGCACTCTTTTGACCAGAGGACTACTTAAATTTGACAG
 ACTGGAGAAGAATGAGGCTGTCTAAAAGAGAAAAAGCCTGATTGTTAAGCAAACTAAAAGCAGAAAAA
 TATAGAAAGATTTACTTAAAGAAATGAAAAAGTTAGATCTCAAGAAATATATAAAGACATTGTGAAG
 AAAAAATTTGTTATGGATATGATTGCCTTTGAAAAAGCCTGTGAAAGACTCAAGATGCTAAAAACAAAGT
 AGCAATTGTGAAAAACAAATTTAGACTTTAAAGTCCCAATGGACTGATAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226040 representing NM_001105659
 Red=Cloning site Green=Tags(s)


MFHGTVEELTSHEEWSHYNENIREGQKDFVFKFNGHLKSMENLQSCISLRVCIFSNNFITDIHPLQS
 CIKLIKLDLHGNQIKSLPNTKFWNGLKNLKLKLLYHDNGFAKLKNICVL SACP TLIAL TMFDCPVSLKKG
 Y RHVLVNSIWPLKALDHHVISDEEIIQNWHLPFRKACNHRLFFNFCPALRKGTYYEEINNKHITSKIN
 AILAHNSPVLIVQRWIRGFLVRKNLSPVFFHKKKQKEKIRGYEAKWIYITKGYEDKLLKDLFFKPETNI
 KGKLAYWKHNIYYVDLKNSSHRKHVSSILCELPKDLGMKSKTSRHLIQKGQSEDEIVDEKLDTSFR
 ISVFKLPIYTSGLKNNAVLREKKQHFFPAYPQPIYTHPKPIIKKDIRLSMKEFFAPQFRAGMKLRTF
 SDIDKYYTEQKKQYHKEKVRVVAQVARERVRVAVNEHLNQKYYATQKLIENKETIQNSLRQVWQNR
 FNYLEKARERKALFLKEKSQKASERLLVQNLNNERTLLTRGLLKIDRLEKNEAVLKEKSLIVKQKLKAEK
 YRKNLLKEMKKVRSQEIYKRHCEEKFVMDMIAFEKACERLQDAKTKVAIVKTNLDFKVPNGLIK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8034_b08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001105659

ORF Size: 1872 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001105659.2](#)


RefSeq ORF: 1875 bp

Locus ID: 127255
 UniProt ID: [A6PVS8](#)
 Cytogenetics: 1p31.1
 MW: 73.5 kDa

Product images:



Circular map for RC226040



Western blot validation of overexpression lysate (Cat# [LY426292]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226040 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).