

Product datasheet for **RC226044**

SHTN1 (NM_001127211) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SHTN1 (NM_001127211) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SHTN1
Synonyms:	KIAA1598; shootin-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC226044 representing NM_001127211
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAACAGCTCGGACGAAGAGAAGCAGCTGCAGCTCATTACCAGTCTGAAGGAGCAAGCAATAGGCGAAT
 ATGAAGACCTTAGAGCAGAGAACCAGAAAAACAAAGGAGAAGTGTGACAAAATTAGGCAAGAACGAGATGA
 AGCCGTTAAAAAACTGGAAGAATTTTCAGAAAAATTTCTCACATGGTCATAGAGGAAGTTAATTTTCATGCAG
 AACCATCTTGAAATAGAGAAGACTTGTCTGAGAAAGTGTGAAGCTTTGGCAACAAAGCTAAATAAAGAAA
 ATAAACGTTGAAAAGAATCAGCATGTTGTACATGGCCAAGCTGGGACCAGATGTAATACTGAAGAGAT
 AAACATTGATGATGAAGATTCGACTACAGACACAGACGGTGCCGCCGAGACTTGTGTCTCAGTACAGTGT
 CAGAAGCAAAATTAAGAACTTCGAGATCAAATTTGTATCTGTTTCTGAGGAGGAAAAGAAGTTTTAGCCATTG
 AGCTGGAAAAATCTCAAGAGCAAACTCGTAGAAGTAATTGAAGAAGTAAATAAAGTTAAACAAGAAAAGAC
 TGTTTTAAATTCAGAAGTTCTTGAACAGAGAAAAGTCTTAGAAAAATGCAATAGAGTGTCCATGTTAGCT
 GTAGAAGAGTATGAGGAGATGCAAGTAAACCTGGAGCTGGAGAAGGACCTTCGAAAGAAAAGCAGAGTCAT
 TTGCACAAGAGATGTTTCATTGAGCAAAAACAGCTAAAGAGACAAAAGCCACCTTCTGCTGCAGAGCTCCAT
 CCCTGATCAGCAGCTTTTGAAAGCTTTAGACGAAAATGCAAACTCACCCAGCACTTGAAGAAGAGAGA
 ATTCAGCATCAACAAAAGGTCAAAGAATTAGAAGAGCACTAGAAAAAGAACTCCACAAAGAAATAC
 ACAACCTCAAACAGCACTGGAGCTTCTAGAGGAAGATAAAAAGGAATTGGAATTGAAATATCAGAAATTC
 TGAAGAGAAAGCCAGAAATTTAAAGCACTCTGTTGATGAACTCCAGAAACGAGTGAACCACTGTGAGAT
 TCAGTACCTCCACCACCTCCTCCTCCACCACCACTTCCCCCTCCACCTCCCAATCCTATCCGATCCCTCA
 TGTCATGATCCGGAACGATCCACCCAGTGGCAGTGGTGCTAAGAAAGAAAAGGCAACTCAACCAGA
 AACAACTGAAGAAGTACAGATCTAAAGAGGCAAGCAGTTGAAGAGATGATGGATAGAATTAAGGGA
 GTTCATCTTAGACCCGTTAATCAGACAGCCAGACCGAAGACAAAGCCAGAATCTTCGAAAGGCTGCGAAA
 GTGCAGTGGATGAACTAAAGGAATACTGGGGACACTTAACAAATCCACTAGTTCAAGAAGCTTAAATC
 CCTTGACCCTGAAAACAGTGAACTGAGTTAGAAAGGATTTTGCCTCGCAGAAAGGTGACAGCAGAAAGCA
 GATAGCAGTAGTCCAAGTGGATATTAGCCACCTCAGAGTCCAAATCCATGCCAGTGTGGGTTCTGTAT
 CCAGTGTAAACAAAACAGCCTTGAACAAGAAAACCTCTGGAGGCAGAAATTAACAGCCCGTCCCCCAAC
 ACCTGAGCCAGGTGAAGGGCCCCGTAATTTGGAAGGATGCACAAGTTCCAAGGTTACGTTTCAGCCTCCC
 AGTAGCATTGGATGCAGGAAAAATACATTGACGGTGAAAAACAAGCCGAACCAAGTTGTAGTTTTAGATC
 CTGTTTCTACACATGAACCCCAACCAAGACCAGGTTGCTGAAAAAGATCCAACCTCAACACAAGGAGGA
 TGAAGGCGAAATTAACAGAAAACAAAGAAGACAGCATTGAAAACGTGAGAGAGACAGACAGCTCCAAC
 TGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226044 representing NM_001127211
 Red=Cloning site Green=Tags(s)


MNSSDEEKQLQLITSLKEQAIGEYEDLRAENQKTKEKCDKIRQERDEAVKKLEEFQKISHMVIEEVNFMQ
 NHLEIEKTCRESAEALATKLNKENKTLKRISMLYMAKLGPDVITEEINIDDEDSTDTDGAATCVSVQC
 QKQIKELRDQIVSVQEEKKILAIENLKSKEVEIEEVNKKVQEKTVLNSEVLEQRKVLEKCNRVSMIA
 VEEYEEMQVNLELEKDLRKKAEFAQEMFIEQNKLRQSHLLQSSIPDQQLKALDENAKLTQQLLEER
 IQHQKVKELEEQLNETLHKEIHNLKQQLLEEDKKELELKYQNSEEKARNLKHSVDELQKRVNQSSE
 SVPPPPPPPPPLPPPPNPIRSLMSMIRKRSHPGSGAKKEKATQPETTEEVTDLKRQAVEEMMDRIKKG
 VHLRPVNQTARPKTKPESSKGCESAVDELKILGTNLKSTSSRLKSLDPENSETELERILRRRKVTAEA
 DSSSPTGILATSESKSMPVLGSVSVTKTALNKKTLAEFNPSPPPTPEPGEPRKLEGCTSSKVTFPQP
 SSIIGCRKKYIDGEKQAEPPVVLDPVSTHEPQTKDQVAEKDPTQHKEDGEIQPENKEDSIENVRETDSSN
 C

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8034_b06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001127211

ORF Size: 1893 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001127211.3</u>
RefSeq ORF:	1896 bp
Locus ID:	57698
UniProt ID:	<u>A0MZ66</u>
Cytogenetics:	10q25.3
MW:	71.5 kDa
Gene Summary:	Involved in the generation of internal asymmetric signals required for neuronal polarization and neurite outgrowth. Mediates netrin-1-induced F-actin-substrate coupling or 'clutch engagement' within the axon growth cone through activation of CDC42, RAC1 and PAK1-dependent signaling pathway, thereby converting the F-actin retrograde flow into traction forces, concomitantly with filopodium extension and axon outgrowth. Plays a role in cytoskeletal organization by regulating the subcellular localization of phosphoinositide 3-kinase (PI3K) activity at the axonal growth cone. Plays also a role in regenerative neurite outgrowth. In the developing cortex, cooperates with KIF20B to promote both the transition from the multipolar to the bipolar stage and the radial migration of cortical neurons from the ventricular zone toward the superficial layer of the neocortex. Involved in the accumulation of phosphatidylinositol 3,4,5-trisphosphate (PIP3) in the growth cone of primary hippocampal neurons.[UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC226044



Western blot validation of overexpression lysate (Cat# [LY426712]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226044 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified SHTN1 protein (Cat# [TP326044]). The protein was produced from HEK293T cells transfected with SHTN1 cDNA clone (Cat# RC226044) using MegaTran 2.0 (Cat# [TT210002]).