

Product datasheet for **RC226049**

c-Myb (MYB) (NM_001130172) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	c-Myb (MYB) (NM_001130172) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	c-Myb
Synonyms:	c-myb; c-myb_CDS; Cmyb; efg
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC226049 representing NM_001130172
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCCCGAAGACCCCGGCACAGCATATATAGCAGTGACGAGGATGATGAGGACTTTGAGATGTGTGACC
ATGACTATGATGGGCTGCTTCCCAAGTCTGAAAAGCGTCACTTGGGGAAAACAAGGTGGACCCGGGAAGA
GGATGAAAAAAGTGAAGAAGCTGGTGAACAGAATGGAACAGATGACTGGAAGTTATTGCCAATTATCTC
CCGAATCGAACAGATGTGCACTGCCAGCACCGATGGCAGAAAAGTAAACCCTGAGCTCATCAAGGGTC
CTTGGACCAAAGAAGAAGATCAGAGAGTGATAGAGCTTGTACAGAAAACGGTCCGAAACGTTGGTCTGT
TATTGCCAAGCACTTAAAGGGGAGAATTGGAACAATGTAGGGAGAGGTGGCATAACCCTTGAATCCA
GAAGTTAAGAAAACCTCCTGGACAGAAGAGGAAGACAGAATTATTTACCAGGCACACAAGAGACTGGGA
ACAGATGGGCAGAAATCGCAAAGCTACTGCCTGGACGAACTGATAATGCTATCAAGAACCCTGGAATTC
TACAATGCGTCGGAAGGTGGAACAGGAAGTTATCTGCAGGAGTCTTCAAAGCCAGCCAGCCAGCAGTG
GCCACAAGCTTCCAGAAGAACAGTCATTTGATGGGTTTTGCTCAGGCTCCGCCTACAGCTCAACTCCCTG
CCTGGCCAGCCACTGTTAACACGACTATTCTATTACCACATTTCTGAAGCACAAAATGTCTCCAG
TCATGTTCCATACCCTGTAGCGTTACATGTAATATAGTCAATGTCCCTCAGCCAGCTGCCCGAGCCATT
CAGAGACACTATAATGATGAAGACCCTGAGAAGGAAAAGCGAATAAAGGAATTAGAATTGCTCCTAATGT
CAACCGAGAATGAGCTAAAAGGACAGCAGACACAGAACCACACATGCAGTACCCCGGTGGCACAGCAC
CACCATTGCCGACCACACCAGACCTCATGGAGACAGTGCACCTGTTTCTGTTGGGAGAACCCACTCC
ACTCCATCTCTGCCAGCGGATCCTGGCTCCCTACCTGAAGAAAGCGCCTCGCCAGCAAGGTGCATGATCG
TCCACCAGGGCACCATTCTGGATAATGTTAAGAACCTTTAGAATTTGCAGAAACACTCCAATTTATAGA
TTCTTTCTTAAACACTTCCAGTAACCATGAAAACCTCAGACTTGGAAATGCCTTCTTAACTTCCACCCC
CTCATTGGTCACAAATTGACTGTTACAACACCATTTTCATAGAGACCAGACTGTGAAAACCTAAAAGGAAA
ATACTGTTTTTAGAACCCAGCTATCAAAGGTCAATCTTAGAAAAGCTCTCAAGAACTCTACACCATT
CAAACATGCAGCTGCAGCTCAAGAAATTAATACGGTCCCCTGAAGATGCTACCTCAGACACCCTCAT
CTAGTAGAAGATCTGCAGGATGTGATCAAACAGGAATCTGATGAATCTGGAATTGTTGCTGAGTTTCAAG
AAAATGGACCACCCTTACTGAAGAAAATCAAACAAGAGGTGGAATCTCCAAGTATAAATCAGGAACTT
CTTCTGCTCACACCCTGGGAAGGGGACAGTCTGAATACCCAAGTGTTCACGCAGACCTCGCCTGTGGCA
GATGCACCGAATATTCTTACAAGCTCCGTTTTAATGGCACCAGCATCAGAAGATGAAGACAATGTTCTCA
AAGCATTTACAGTACCTAAAACAGTCCCTGGCGAGCCCTTGCAGCCTTGTAGCAGTACCTGGGAACC
TGATCCTGTGGAAGATGGAGGAGCAGATGACATCTTCCAGTCAAGCTCGTAAATACGTGAATGCATTC
TCAGCCCGGACGCTGGTCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226049 representing NM_001130172
Red=Cloning site Green=Tags(s)


```
MARRPRHSIYSSDEDDDFEMCDHDYDGLLPKSGKRHLGKTRWTREEDEKLKLVQNGTDDWKVIANYL  
PNRTDVQCQHRWQKVLNPELIKGPWTKEEDQRVIELVQKYGPKRWSVIAKHLKGRIGKQCRERWHNHLNP  
EVKKTSWTEEDRIIYQAHKRLGNRWAEIAKLLPGRTDNAIKNHWNSTMRKVEQEGYLQESSKASQPAV  
ATSFQKNSHLMGFAPPTAQLPATGQPTVNDYSYHISEAQNVSSHVPYPVALHVNIIVNPQAAAAI  
QRHYNDEDPEKEKRIKELELLLMSTENELKGQQTQNHSTCSYPGWSTTIADHTRPHGDSAPVSLGEHHS  
TPSLPADPGSLPEESASPARCMIHQGTILDNVKNLLEFAETLQFIDSFLNTSSNHENSLEMPSTSTP  
LIGHKLTVTTPFHRDQTVKTQKENTVFRTPAIKRSILESSPRTPFPFKHALAAQEKYGLKMLPQTPSH  
LVEDLQDVIKQESDESGIVAEFQENGPPLKKIKQEVESPTDKSGNFFCSHHWEGDSLNTQLFTQTSPVA  
DAPNILTSSVLMAPASEDEDNVLAFTVPKNRSASPLQPCSSTWEPASCGKMEEQMTSSSQARKYVNAF  
SARTLVM
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8070_g11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001130172

ORF Size: 1911 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130172.2](#)

RefSeq ORF: 1914 bp

Locus ID: 4602

UniProt ID: [P10242](#)


Cytogenetics: 6q23.3

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Stem cell - Pluripotency, Transcription Factors

MW: 71.9 kDa

Gene Summary: This gene encodes a protein with three HTH DNA-binding domains that functions as a transcription regulator. This protein plays an essential role in the regulation of hematopoiesis. This gene may be aberrantly expressed or rearranged or undergo translocation in leukemias and lymphomas, and is considered to be an oncogene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]

Product images:



Circular map for RC226049