

## Product datasheet for **RC226061**

### **PTPN12 (NM\_001131009) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PTPN12 (NM_001131009) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PTPN12
Synonyms:	PTP-PEST; PTPG1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC226061 representing NM\_001131009  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGCCTGCCGAGAATTTGAGATGGGAAGGAAAAATGTGAGCGCTATTGGCCTTTGTATGGAGAAGACC  
CCATAACGTTTGCACCATTTAAAATTTCTTGTGAGGATGAACAAGCAAGAACAGACTACTTCATCAGGAC  
ACTCTTACTTGAATTTCAAATGAATCTCGTAGGCTGTATCAGTTTCATTATGTGAAGTGGCCAGACCAT  
GATGTTCCCTTCATCATTGATTCTATTCTGGACATGATAAGCTTAATGAGGAAATATCAAGAACATGAAG  
ATGTTCCATTTGATTTCATTGCAGTGCAGGCTGTGGAAGAACAGGTGCCATTTGTGCCATAGATTATAC  
GTGGAATTTACTAAAAGCTGGGAAAATACCAGAGGAATTAATGTATTTAATTAACAAGAAATGAGA  
ACACAAAGGCATTCTGCAGTACAAACAAAGGAGCAATATGAACTTGTTTCATAGAGCTATTGCCAACTGT  
TTGAAAAACAGCTACAACATATGAAATTCATGGAGCTCAGAAAATGCTGATGGAGTGAATGAAATTA  
CACTGAAAACATGGTCACTCCATAGAGCTGAAAACAAGATTCTCCTCCTCCAAAACCAAGGACC  
CGCAGTTGCCCTTGTGAAGGGGATGCTAAAGAAGAAACTGCAGCCACCGGAACCTCATCAGTGCCAC  
CCATCTTGACACCTTCTCCCCCTCAGCTTTTCCAACAGTCACTACTGTGTGGCAGGACAATGATAGATA  
CCATCCAAAGCCAGTGTTCATATGTTTCATCAGAACAACATTCAGCAGACCTCAACAGAACTATAGT  
AAATCAACAGAATTCAGGGAAAAATGAATCAACAATGAACAGATAGATAAAAAATGGAAACGAAAT  
TAAGTTTTGAGATTAAGAAGTCCCTCTCCAAGAGGGACCAAAAAGTTTTGATGGGAACACACTTTTGA  
TAGGGGACATGCAATTAATAAATCTGCTTACCTTGATAGCTGATAAAATCTCTAAGCCACAGGAA  
TTAAGTTCAGATCTAAATGTCGGTGATACTTCCAGAATTCTTGTGTGGACTGCAGTGAACACAATCAA  
ACAAAGTTTCAGTTACTCCACCAGAAGAATCCAGAATTCAGACACACCTCCAAGGCCAGCCGCTTGC  
TCTTGATGAGAAAAGGACATGTAACGTGGTCATTTTCATGGACCTGAAAATGCCATACCCATACCTGATTTA  
TCTGAAGCAATTCCTCAGATATCAACTATCAAACTAGGAAAATGTGAGTTTAAACCAAGTCTACAA  
CACAAGTTGAAACACCTGATCTTGTGGATCATGATAAACAATTCACCACTCTTCAGAACACCCCTCAGTTT  
TACTAATCCACTTCACTCTGATGACTCAGACTCAGATGAAAGAACTCTGATGGTGTGTGACCCAGAAT  
AAAATAATTTCAACAGCAAGTCCACAGTTTCTGCTGCCACTAGTACTGAAAGCATTCTACTAGGA  
AAGTATTGCCAATGTCCATTGCTAGACATAATATAGCAGGAACAACACATTCAGGTGCTGAAAAAGATGT  
TGATGTTAGTGAAGATTCACCTCCTCCCCTACCTGAAAGAACTCCTGAATCGTTTGTGTTAGCAAGTGAA  
CATAATACACCTGTAAGATCGGAATGGAGTGAACCTCAAAGTCAAGAACGATCTGAACAAAAAAGTCTG  
AAGGCTTGATAACCTCTGAAAATGAGAAATGTGATCATCCAGCGGAGGTATTCATATGAAATGTGCAT  
AGAATGTTCCACTACTTTCAAGTACAAGAGAGAAACAATATCAGAAAATCCAACAGAAGCCACAGATATT  
GGTTTTGGTAATCGATGTGGAACCCAAAGGACCAAGAGATCCACCTTCAGAATGGACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226061 representing NM\_001131009  
 Red=Cloning site Green=Tags(s)


MACREFEMGRKKCERYWPLYGEDPITFAPFKISCEDEQARTDYFIRTLLEFQNESRRLYQFHVYNWPDH  
 DVPSSFDISLDMISLMRKYQEHEDVPICIHCSAGCGRTGAICAIDYTWLLKAGIPEEFNVFNLIQEMR  
 TQRHSAVQTKEQYELVHRAIAQLFEKQLQLYEIHGAQKIADGVNEINTENMVSSIEPEKQDSPPKPPRT  
 RSCLVEGDAKEEILQPPEPHVPPILTPSPPSAFPTVTTVWQDNDRYHPKPVLMHVSSEQHSADLNRNYS  
 KSTELPGKNESTIEQIDKKLERNLSFEIKKVPLQEGPKSFDGNTLLNRGHAIKIKSASPCIADKISKPQE  
 LSSDLNVGDSQNSCVDCSVTQSNKVSVPPEESQNSDTPPRPDRLPLDEKGHVTSFHGPNNAIPIDL  
 SEGSSDINYQTRKTVSLTPSPTTQVETPDLVDHDNTSPLFRTPLSFTNPLHSDDSDSDERNSDGAVTQN  
 KTNISTASATVSAATSTESISTRKVLPMIARHNIAGTTHSGAEKDVDVSEDSPPPLPERTPESEVLASE  
 HNTPVRSWSELQSQRSEQKKSEGLITSENEKCDHPAGGIHYEMCIECPPTFSDKREQISENPTEATDI  
 GFGNRCGKPKGPRDPPSEWT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001131009

ORF Size: 1950 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001131009.2](#)

RefSeq Size:	3026 bp
RefSeq ORF:	1953 bp
Locus ID:	5782
UniProt ID:	<a href="#">Q05209</a>
Cytogenetics:	7q11.23
Protein Families:	Druggable Genome, Phosphatase
MW:	72.8 kDa

**Gene Summary:** The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains a C-terminal PEST motif, which serves as a protein-protein interaction domain, and may regulate protein intracellular half-life. This PTP was found to bind and dephosphorylate the product of the oncogene c-ABL and thus may play a role in oncogenesis. This PTP was also shown to interact with, and dephosphorylate, various products related to cytoskeletal structure and cell adhesion, such as p130 (Cas), CAKbeta/PTK2B, PSTPIP1, and paxillin. This suggests it has a regulatory role in controlling cell shape and mobility. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008]