

## Product datasheet for **RC226075L1V**

### **LRR32 (NM\_001128922) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	LRR32 (NM_001128922) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LRR32
Synonyms:	CPPRDD; D11S833E; GARP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001128922
ORF Size:	1986 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226075).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001128922.1</a> , <a href="#">NP_001122394.1</a>
RefSeq Size:	4222 bp
RefSeq ORF:	1989 bp
Locus ID:	2615
UniProt ID:	<a href="#">Q14392</a>
Cytogenetics:	11q13.5
Protein Families:	Transmembrane
MW:	72 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a type I membrane protein which contains 20 leucine-rich repeats. Alterations in the chromosomal region 11q13-11q14 are involved in several pathologies. [provided by RefSeq, Jul 2008]