

Product datasheet for **RC226079**

NEDD1 (NM_001135175) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NEDD1 (NM_001135175) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NEDD1
Synonyms:	GCP-WD; TUBGCP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC226079 representing NM_001135175
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCATTTTACAGGCGCAGTCATGCAGGAAAACCTCAGATTTGCTTCATCAGGAGATGATATTTAAATAT
GGGATGCTTCATCTATGACATTGGTGGATAAATTCAACCCACACACATCACCACATGGAATCAGCTCAAT
ATGTTGGAGCAGCAATAATAACTTTTTAGTAACAGCATCTCCAGTGGCGACAAAATAGTTGTCTCAAGT
TGCAAATGTAACCTGTTCCACTTTAGAGCTTGTGAAGGGCAAAGCAGACATGTGTCAATTTAAATT
CTACATCTATGTATTTGGTAAGCGGAGGCCATAAACAACCTGTTAATATTTGGGATTTAAATCAAAAAG
AGTTCATCGATCTCTAAGGATCATAAAGATCAAGTAACCTGTGTAAACATAACAATTGGAATGATTGCTAC
ATTGCTTCTGGATCTCTAGTGGTGAATTTATTTACACAGTGAACCACTAATTTATCTAGTACTCCTT
TTGGCCATGGTAGTAACCACTGTTCCGCACTGAAGTACTCCTGTTAAGAAATCACTACTGGCGAG
TGTTTCGGATAATGGAATAGTAACCTCTCGGGATGTAATAGTCAGAGTCCATACCATAACTTTGACAGT
GTACACAAAGCTCCAGCGTCAGGCATCTGTTTTCTCCTGTCAATGAATTGCTCTTTGTAACCATAGGCT
TGGATAAAAAGAACATCCTCTATGACACTTCAAGTAAGAAGCTAGTGAAGAACTTTAGTGGCTGACTCC
TCTAAGTCCGGTAGATTTTATGACCTGATGGAGCCACTTTGGCTATTGGATCTTCCCGGGGAAAAATAT
CAATATGATTTAAGAATGTTGAAATCACCAGTTAAGACCATCAGTGCTCACAAGACATCTGTGCAGTGTA
TAGCATTTCACTACTCCACTGTTCTTACTAAGTCAAGTTAAATAAAGGCTGTTCAAATAAGCCACAAC
AGTGAACAAACGAAGTGTAAATGTGAATGCTGCTAGTGGAGGAGTTCAGAATCCGGAATGTGCAGAGAA
GCACCTGCCACGTCATTGCCACAGTTCTACCACAACCTATGACATCAGCTATGGGAAAAGGAACAGTTG
CTGTTCAAGAAAAGCAGGTTTGCCTCGAAGCATAAACACAGACACTTTATCTAAGGAAAACAGACAGTGG
AAAAAATCAGGATTTCTCCAGCTTTGATGATACTGGGAAAAGTAGTTTAGGTGACATGTTCTCACCTATC
AGAGATGATGCTGTAGTTAACAAGGGAAGTGTAGTCCATAGGCAAAGGAGATGGCTTTGACTTTCTAC
CGCAGTTGAACCTCAGTGTTCCTCCAAGAAAAATCCAGTAACTTCAAGTACTTCAAGTATTGCAATCTAG
TCCTCTTAATGTTTTTATGGGATCTCCAGGAAAAGAGGAAAATGAAAACCGTGATCTAACAGCTGAGTCT
AAGAAAATATATGGGAAAACAGGAATCTAAAGACTCCTTCAAACAGTTAGCAAAGTTGGTCACATCTG
GTGCTGAAAGTGGAAATCTAAATACCTCTCCATCATCTAACCAACAAGAAATCTGAGAAATTTGAAAA
GCCAGAGAATGAAATGAAGCCAGTTGATATGTGAACCCCAATCAATGGATCCTCAACTCCAAATCCA
AAGATAGCATCTTCTGCTACTGCTGGAGTTGCCAGTTCCTCAGAAAAATAGCCGACAGCATTGGAA
ATAACCGCAAAATGCACCATGACTTCCATTCAAATTCGTTTTATTAGAACATGATACAGGAAACGTT
GGATGACTTTAGAGAAGCATGCCATAGGACATTGTGAATTTGCAAGTGGAGATGATTAACAGTTTCAT
ATGCAACTGAATGAAATGCATTCTTTGCTGAAAAGATACTCAGTGAATGAAGTTTAGTGGCTGAAATG
AAAGACTACGAGAAGAAAACAAAAGATTACGGGCCCACTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226079 representing NM_001135175
Red=Cloning site Green=Tags(s)

```
MHFTGAVMQENLRFASSGDDIKIWDASSMTLVDFKPNHTSPHGISSICWSSNNFLVTASSSGDKIVVSS
CKCKPVPLLELAEGQKQTCVNLNSTSMYLVSGGLNNTVNIWDLKSKRVHRSLKDHKDQVTCVTYNWDCY
IASGSLSGEIIILHSVTTLNSTPFGHGSNQSVRHLKYSLFKKSLLGSVSDNGIVTLWDVNSQSPYHNFDS
VHKAPASGICFSPVNEELLFVTIGLDKRIILYDTSKKLVKTLVADTPLTAVDFMPDGATLAIGSSRGKIY
QYDLRMLKSPVKTISAHKTSVQCI AFQYSTVLTKSSLNKGCSNKPTTVNKRSVNVNAASGGVQNSGIVRE
APATSIATVLPQPMTSAMGKGTVAVQEKAAGLPRINTDTLSKETDSGKNQDFSSFDGKSSLDGMFSPI
RDDAVVNKGSDESIGKGDGDFLPLQLNSVFPPRKNPVTSSVSLHSSPLNVFMGSPGKEENENRDLTAES
KKIYMGKQESKDSFKQLAKLVTSAGESGNLNTSPSSNQTRNSEKFEKPEIEAQLICEPPINGSSTPNP
KIASSVTAGVASSLSEKIADSIGNNRQNAPLTSIQIRFIQNIQETLDDFREACHRDIVNLQVEMIKQFH
MQLNEMHSLLEYSVNEGLVAEIERLREENKRLRAHF
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001135175

ORF Size: 2001 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135175.1](#), [NP_001128647.1](#)

RefSeq Size: 3805 bp

RefSeq ORF: 2004 bp

Locus ID: 121441


UniProt ID: [Q8NHV4](#)

Cytogenetics: 12q23.1

MW: 72.7 kDa

Gene Summary: Required for mitosis progression. Promotes the nucleation of microtubules from the spindle. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC226079