

Product datasheet for **RC226113**

PSD95 (DLG4) (NM_001128827) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PSD95 (DLG4) (NM_001128827) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PSD95
Synonyms:	MRD62; PSD95; SAP-90; SAP90
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC226113 representing NM_001128827
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGACTGTCTCTGTATAGTGACAACCAAGAAATACCGCTACCAAGATGAAGACACGCCCCCTCTGGAGC
 ACAGCCCGGCCACCTCCCAACCAGGCCAATTCTCCCCAGTGATTGTCAACACAGATACCCTAGAAGC
 CCCAGGATATGTGAACGGGACCGAGGGGAGATGGAATACGAGGAAATCACATTGAAAAGGGTAACTCA
 GGTCTGGGCTTCAGCATCGCAGGTGGCACTGACAACCCACACATCGGTGACGACCCATCCATTTTCATCA
 CCAAGATCATTCTGGTGGGCTGCGGCCAGGATGGCCGCTCAGGGTCAACGACAGCATCCTGTTTGT
 AAATGAAGTGGACGTGCGCGAGGTGACCCACTCAGCGCGGTGGAAGCCCTCAAAGAGGCAGGCTCCATC
 GTTCGCTCTATGTCATGCGCCGAAGCCCCGGCTGAGAAGTTCATGGAGATCAAGTTCATCAAGGGC
 CTAAGGTCTTGGCTTCAGCATCGCAGGGGCGTAGGGAACCAGCACATCCAGGAGATAATAGCATCTA
 TGTAACAAAGATCATCGAAGGGGTGCTGCCACAAGGATGGGAGTTGCAGATTGGAGACAAGATCCTG
 GCGGTCAACAGTGTGGGCTAGAGGACGTCATGCATGAAGATGCTGTGGCAGCCCTGAAGAACACGTATG
 ATGTTGTCTACCTAAAGTGGCCAAGCCAGCAATGCCTACCTGAGTGACAGCTATGCTCCCCCAGACAT
 CACAACCTCTTATTTCCAGCACCTGGACAATGAGATCAGTCACAGCAGCTACCTGGGCACCGACTACCC
 ACAGCCATGACCCCACTTCCCCGCGCTACTCTCCAGTGGCCAAGGACCTGCTCGGGGAGGAAGACA
 TCCCCGAGAACCAGGGCAATTGTGATCCACCGGGCTCCACGGGCTGGGCTTCAACATCGTGGGTGG
 CGAGGACGGTGAAGGCATCTTCTCTCTTTATCTGGCCGGGGCCCTGCAGACCTCAGTGGGGAGCTG
 CGAAGGGGGACAGATCCTGTCGGTCAACGGTGTGGACCTCCGAAATGCCAGCCATGAGCAGGCTGCCA
 TTGCCCTGAAGAATGCGGGTCAGACGGTCAGATCATCGCTCAGTATAAACCAGAAGAGTACAGCCGATT
 CGAGGCCAAGATCCACGACCTTCGGGAACAGTCAATGAACAGCAGCCTGGGCTCAGGACTGCGTCCCTG
 CGGAGCAACCCAAAAGGGTTTCTACATCAGGGCCCTGTTTGATTACGACAAGACCAAGGACTGCGGCT
 TCCTGAGCCAGGCCCTGAGCTTCCGCTTTGGGATGTGCTGCATGTCATCGATGCTAGTGTGAGGAGTG
 GTGGCAGGCACGGCGGTCCACTCTGACAGTGAGACCGACGACATTGGGTTTATCCCCAGCAAACGGCGG
 GTTGAGCGACGAGAGTGGTCAAGTTAAAGGCCAAGGACTGGGGCTCCAGCTCTGGATCGCAGGGTCGAG
 AAGACTCGTTCTGAGCTACGAGACAGTGACGCAGATGGAAGTGCATATGCTCGCCCCATCATCATCT
 TGGGCCACCAAGGACCGCCCAACGATGATCTTCTCTCCGAGTCCCCGACAAGTTGGATCCTGTGTT
 CCCCATACGACACGGCCCAAGCGGGAGTATGAGATAGATGGCCGGGATTACCACTTTGTGTGTCGCCGG
 AGAAAATGGAGAAGGACATTAGGCGCACAAAGTTCATTGAGGCCGCCAGTACAACAGCCACCTCTATGG
 GACCAGCTCCAGTCCGTGCGAGAGGTGGCAGAGCAGGGGAAGCACTGCATCCTCGATGTCTCGGCAAT
 GCCGTGCGGGCGCTGCAGGCGGCCACCTGCACCCCATCGCCATCTTATCCGCCCCCGCTCCCTGGAGA
 ATGTGCTAGAGATTAACAAGCGGATCACAGAGGAGCAAGCCCGCAAAGCCTTCGACAGAGCCACCAAGCT
 GGAGCAGGAGTTCACAGAGTCTTCTCAGCCATCGTGGAGGGTGACAGCTTTGAGGAGATCTACCACAAG
 GTGAAGCGTGTATCGAGGACCTCTCAGGCCCTACATCTGGGTTCCAGCCCGAGAGAGACTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226113 representing NM_001128827
 Red=Cloning site Green=Tags(s)


MDCLCIVTTKKYRYQDEDTPPLEHSPAHLPNQANSPPVIVNTDTLEAPGYVNGTEGEMEEYEEITLERNGS
 GLGFSIAGGTDNPHIGDDPSIFITKIIPGGAAAQDGLRVNDSILFVNEVDVREVTHSAAVEALKEAGSI
 VRLYVMRRKPPAEKVMEIKLIKPKGLGF SIAGGVGNQHIPGDNSIYVTKIEGGAHKDGRQLIGDKIL
 AVNSVGLLEDVMHEDAVAALKNTYDVVYLKVAKPSNAYLSDSYAPPDITTSYSQHLDNEISHSSYLGTDYP
 TAMTPTSPRRYSPVAKDLLGEEDIPREPRRIVIHGSGTGLGFNI VGGEDGEGIFISFILAGGPADLSGEL
 RKGDQILSVNGVDLRNASHEQAAIALKNAGQVTIIAQYKPEEYSRFEAKIHDLREQLMNSSLSGSGTASL
 RSNPKRGFYIRALFDYDKTKDCGFLSQALSFRFGDVLHVIDASDEEWQARRVHSDSETDDIGFIPSKRR
 VERREWSRLKAKDWGSSSGSQGREDSVLSYETVTQMEVHYARPIIILGPTKDRANDDLLSEFPDKFGSCV
 PHTTRPKREYIDGRDYHFVSSREKMEKDIQAHKFI EAGQYNSHLYGTSVQSVREVAEQGKHCILDV SAN
 AVRRLQAAHLHPAIFIRPRSLNVLEINKRITTEEQARKAFDRATKLEQEFTECFSAIVEGDSFEEIYHK
 VKRVIEDLSGPYIWWPARERL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001128827

ORF Size: 2163 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128827.3](#)

RefSeq ORF: 2166 bp

Locus ID: 1742

UniProt ID: [P78352](#)

Cytogenetics: 17p13.1


Protein Families: Druggable Genome

Protein Pathways: Huntington's disease

MW: 79.9 kDa

Gene Summary: This gene encodes a member of the membrane-associated guanylate kinase (MAGUK) family. It heteromultimerizes with another MAGUK protein, DLG2, and is recruited into NMDA receptor and potassium channel clusters. These two MAGUK proteins may interact at postsynaptic sites to form a multimeric scaffold for the clustering of receptors, ion channels, and associated signaling proteins. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC226113