

Product datasheet for **RC226131L3V**

ISLR2 (NM_001130137) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ISLR2 (NM_001130137) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ISLR2
Synonyms:	LINX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130137
ORF Size:	2235 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226131).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130137.1 , NP_001123609.1
RefSeq Size:	4569 bp
RefSeq ORF:	2238 bp
Locus ID:	57611
UniProt ID:	Q6UXK2
Cytogenetics:	15q24.1
Protein Families:	Druggable Genome, Transmembrane
MW:	79 kDa



[View online »](#)

Gene Summary: Required for axon extension during neural development.[UniProtKB/Swiss-Prot Function]