

Product datasheet for **RC226135**

CERT1 (NM_001130105) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CERT1 (NM_001130105) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CERT1
Synonyms:	CERT; CERTL; COL4A3BP; GPBP; MRD34; STARD11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC226135 representing NM_001130105
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGCACAGCTGCATCCCTACCCCGCCTCTCCTTTCTCCGCTCCTCTGCTTTTCTACCCGTCGTCA
 CCCGGGAGAGCCGGAGGGGGCTAAGTTCGGGTGGCAGCGCCGGGCGCAACGCAGGGGTACGGCGACGGC
 GGCGGGCGCTGACGGCTGGAAGGGTAGGCTTCCTTACCGCTCGTCTCCTTCTCGCTCCGCTCGGTGT
 CAGGCGCGGGCGGCGCGGGGGGACTTCGTCCCTCCTGCTCCCCCACACGGAGCGGGCAC
 TCTTCGTTCCGCATCCCCGACCTTACCCGAGGACTGGGCGCTCCTCCGGCGCAGCTGAGGGAGC
 GGGGGCCGGTCTCCTGCTCGTTGTGAGCCTCCATGTCGATAATCAGAGCTGGAACCTCGTCGGGCTCG
 GAGGAGGATCCAGAGACGGAGTCTGGGCGCCTGTGGAGCGCTCGGGGTCTCAGTAAGTGACAAACT
 ACATTCATGGGTGGCAGGATCGTTGGTAGTTTTGAAAAATAATGCTCTGAGTTACTACAACTGAAGA
 TGAACAGAGTATGGCTGCAGAGGATCCATCTGTCTTAGCAAGGCTGTATCACACCTCAGATTTTGAT
 GAATGTCGATTTGATATTAGTGTAAATGATAGTGTGGTATCTTCGTGCTCAGGATCCAGATCATAGAC
 AGCAATGGATAGATGCCATTGAACAGCAAGACTGAATCTGGATATGGATCTGAATCCAGCTTGCCTCG
 ACATGGCTCAATGGTGTCCCTGGTGTCTGGAGCAAGTGGCTACTCTGCAACATCCACCTTTCATTCAAG
 AAAGGCCACAGTTTACGTGAGAAGTTGGCTGAAATGAAAACATTTAGAGACATCTTATGTAGACAAGTTG
 ACACGCTACAGAAGTACTTTGATGCCTGTGCTGATGCTGTCTTAAGGATGAACCTCAAAGGGATAAAGT
 GGTAGAAGATGATGAAGTACTTTTCTACAACCGTCTGATGGTACTTCTTCATAGTACCAACGGC
 AATAAAGAAAAGTTATTTCCACATGTGACACAAAAGGAATTAATGGTATAGACTTTAAAGGGAAAGCGA
 TAACTTTTAAAGCAACTACTGCTGGAATCCTTGCAACACTTTCTCATTGTATTGAACTAATGGTTAAACG
 TGAGGACAGCTGGCAGAAGAGACTGGATAAGGAAACTGAGAAGAAAAGAAGAACAGAGGAAGCATATAAA
 AATGCAATGACAGAACTTAAGAAAAATCCCACTTTGGAGGACCAGATTATGAAGAAGGCCCTAACAGTC
 TGATTAATGAAGAAGAGTTCTTTGATGCTGTTGAAGCTGCTCTTGACAGACAAGATAAAATAGAAGAACA
 GTCACAGAGTGAAAAGGTGAGATTACATTGGCCTACATCCTTGGCCTCTGGAGATGCCTTTTCTCTGTG
 GGGACACATAGATTTGTCCAAAAGCCCTATAGTCGCTCTTCTCCATGTCTTCCATTGATCTAGTCAGTG
 CCTCTGATGATGTTACAGATTCAGCTCCAGGTTGAAGAGATGGTGCAGAACCACATGACTTACTCATT
 ACAGGATGTAGGCGGAGATGCCAATTGGCAGTTGGTTGTAGAAGAAGGAGAAATGAAGGTATACAGAAGA
 GAAGTAGAAGAAAATGGGATTGTTCTGGATCCTTTAAAAGCTACCCATGCAGTTAAAGGCGTCACAGGAC
 ATGAAGTCTGCAATTATTTCTGGAATGTTGACGTTTCGCAATGACTGGGAAACAACATATAGAAAACCTTCA
 TGTGGTGGAAACATTAGCTGATAATGCAATCATATTTATCAAACACACAAGAGGGTGTGGCCTGCTTCT
 CAGCGAGACGTATTATATCTTTCTGTCATTGAAAAGATACCAGCCTTGACTGAAAATGACCCGAAAACCT
 GGATAGTTTGTAAATTTTCTGTGGATCATGACAGTGTCTCTAAACAACCGATGTGTCCGTGCCAAAAT
 AAATGTTGCTATGATTTGTCAAACCTTGGTAAGCCACAGAGGGAAACCAGGAAATTAGCAGGGACAAC
 ATTCTATGCAAGATTACATATGTAGCTAATGTGAACCTGGAGGATGGGCACCAGCCTCAGTGTTAAGGG
 CAGTGGCAAAGCGAGAGTATCCTAAATTTCTAAAACGTTTTACTTCTTACGTCCAAGAAAAAAGTGCAGG
 AAAGCCTATTTTGTTT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC226135 representing NM_001130105
Red=Cloning site Green=Tags(s)

MQHSCIPTPPSPFSAPPAFLPVVTRESRRGLSSGGSAGRAGVTATAAADGWKGRPLPSPLVLLPRSARC
 QARRRRGGRTSSL LLLPPTPERALFASPPDPSPRGLGASSGAAEGAGAGLLGCRASMSDNQSWNSSGS
 EEDPETESGPPVERCGVL SKWTNYIHGWQDRWVVLKNNALSYKSEDETEYGCRGSICLSKAVITPHDFD
 ECRFDISVNDVWYLRAQDPDHRQQWIDAIEQHKTESGYGSESSLRRHGSMVSLVSGASGYSATSTSSFK
 KGHSLREKLAEMETFRDILCRQVDTLQKYFDACADAVSKDELQRDKVVEDEDDDFPTTRSDGDFLHSTNG
 NKEKLFPHVTPKINGIDFKGEAITFKATTAGILATLSHCIELMKREDSWQKRLDKETEKRRTEEAYK
 NAMTELKKKSHFGGPDYEEGPNSLINEEEFFDAVEAALDRQDKIEEQSQSEKVRHLHWPTSLPSGDAFSSV
 GTHRFVQKPYSRSSMSIDLVSASDDVHRFSSQVEEMVQNHMTYSLQDVGGDANWQLVVEEGEMKVYRR
 EVEENGI VLDPLKATHAVKGV TGHEVCNYFWNVDRNDWETT IENFHVETLADNAII IYQTHKRVWPAS
 QRDVLYLSVIRKIPAL TENDPETWIVCNF SVDHDSAPLNNRCVRAKINVAMICQTLVSPPEGNQEI SRDN
 ILCKITYVANVNP GGWAPASVLR AVAKREYKFLKRF TSYVQEKTAGKPI LF

SGPTRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001130105

ORF Size: 2256 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130105.1](#), [NP_001123577.1](#)

RefSeq ORF: 2259 bp

Locus ID: 10087


UniProt ID: [Q9Y5P4](#)

Cytogenetics: 5q13.3

MW: 83.5 kDa

Gene Summary: This gene encodes a kinase that specifically phosphorylates the N-terminal region of the non-collagenous domain of the alpha 3 chain of type IV collagen, known as the Goodpasture antigen. Goodpasture disease is the result of an autoimmune response directed at this antigen. One isoform of this protein is also involved in ceramide intracellular transport. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC226135