

## Product datasheet for **RC226135L3V**

### **COL4A3BP (CERT1) (NM\_001130105) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	COL4A3BP (CERT1) (NM_001130105) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CERT1
Synonyms:	CERT; CERTL; COL4A3BP; GPBP; MRD34; STARD11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130105
ORF Size:	2256 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226135).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001130105.1</a> , <a href="#">NP_001123577.1</a>
RefSeq ORF:	2259 bp
Locus ID:	10087
UniProt ID:	<a href="#">Q9Y5P4</a>
Cytogenetics:	5q13.3
MW:	83.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a kinase that specifically phosphorylates the N-terminal region of the non-collagenous domain of the alpha 3 chain of type IV collagen, known as the Goodpasture antigen. Goodpasture disease is the result of an autoimmune response directed at this antigen. One isoform of this protein is also involved in ceramide intracellular transport. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]