

## Product datasheet for **RC226141L2V**

### PWWP2A (NM\_001130864) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PWWP2A (NM_001130864) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PWWP2A
Synonyms:	MST101
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001130864
ORF Size:	2265 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226141).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001130864.1</a>
RefSeq ORF:	2268 bp
Locus ID:	114825
UniProt ID:	<a href="#">Q96N64</a>
Cytogenetics:	5q33.3
MW:	81.8 kDa
Gene Summary:	H2A.Z-specific chromatin binding protein which may play an important role in the neural crest stem cell migration and differentiation during early development. Also required for proper mitosis progression.[UniProtKB/Swiss-Prot Function]


[View online »](#)