

Product datasheet for **RC226145L4V**

STRIP2 (NM_001134336) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	STRIP2 (NM_001134336) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STRIP2
Synonyms:	FAM40B; FAR11B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001134336
ORF Size:	2274 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226145).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001134336.1
RefSeq ORF:	2277 bp
Locus ID:	57464
UniProt ID:	Q9ULQ0
Cytogenetics:	7q32.1
MW:	85.8 kDa
Gene Summary:	Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape.[UniProtKB/Swiss-Prot Function]


[View online »](#)