

Product datasheet for **RC226154L3V**

SH3D19 (NM_001128923) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SH3D19 (NM_001128923) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SH3D19
Synonyms:	EBP; Eve-1; EVE1; Kryn; SH3P19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001128923
ORF Size:	2301 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226154).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001128923.1
RefSeq Size:	5210 bp
RefSeq ORF:	2304 bp
Locus ID:	152503
UniProt ID:	Q5HYK7
Cytogenetics:	4q31.3
MW:	84.1 kDa


[View online »](#)

Gene Summary:

This gene encodes a multiple SH3 domain-containing protein, which interacts with other proteins, such as EBP and members of ADAM family, via the SH3 domains. This protein may be involved in suppression of Ras-induced cellular transformation and Ras-mediated activation of ELK1 by EBP, and regulation of ADAM proteins in the signaling of EGFR-ligand shedding. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]