

## Product datasheet for **RC226162L1V**

### **cPLA2 beta (PLA2G4B) (NM\_001114633) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	cPLA2 beta (PLA2G4B) (NM_001114633) Human Tagged ORF Clone Lentiviral Particle
Symbol:	cPLA2 beta
Synonyms:	cPLA2-beta; HsT16992
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001114633
ORF Size:	2343 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226162).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001114633.1</a> , <a href="#">NP_001108105.1</a>
RefSeq ORF:	2346 bp
Locus ID:	100137049
UniProt ID:	<a href="#">P0C869</a>
Cytogenetics:	15q15.1
MW:	87.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the cytosolic phospholipase A2 protein family. Phospholipase A2 enzymes hydrolyze the sn-2 bond of phospholipids, releasing lysophospholipids and fatty acids. This enzyme may be associated with mitochondria and early endosomes. Most tissues also express read-through transcripts from the upstream gene into this gene, some of which may encode fusion proteins combining the N-terminus of the upstream gene including its JmjC domain with the almost complete coding region of this gene, including the C2 and cytoplasmic phospholipase A2 domains. [provided by RefSeq, Jul 2008]