

## Product datasheet for **RC226164L2V**

### HHIPL1 (NM\_001127258) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HHIPL1 (NM_001127258) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HHIPL1
Synonyms:	KIAA1822; UNQ9245
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001127258
ORF Size:	2346 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226164).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001127258.1</a>
RefSeq ORF:	2349 bp
Locus ID:	84439
UniProt ID:	<a href="#">Q96JK4</a>
Cytogenetics:	14q32.2
Protein Families:	Druggable Genome
MW:	86.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that belongs to the glucose/sorbosone dehydrogenase family. The encoded protein also contains a domain that binds folate and reduced folic acid derivatives. [provided by RefSeq, Jul 2016]