

Product datasheet for **RC226170**

KIF9 (NM_001134878) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIF9 (NM_001134878) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KIF9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC226170 representing NM_001134878
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGTACTAGGAAAAAGTTCATGCATTTGTCCGTGTCAAACCCACCGATGACTTTGCTCATGAAATGA
 TCAGATACGGAGATGACAAAAGAAGCATTGATATTCACCTTAAAAAAGACATTCGGAGAGGAGTTGTCAA
 TAACCAACAGACAGACTGGTCGTTTAAGTTGGATGGAGTTCTTCACGATGCCTCCAGGACTTGGTTTAT
 GAGACAGTTGCAAAGGATGTGGTTTCTCAGGCCCTCGATGGCTATAATGGCACCATCATGTGTTATGGGC
 AGACGGGAGCTGGCAAGACATACACCATGATGGGGCAACTGAGAATTACAAGCACCGGGGGATCTCCC
 TCGTGCCCTGCAGCAGGTTTTTAGGATGATCGAAGAACGCCCCACACATGCCATCACTGTGCGTGTTC
 TACTTGAAATCTATAATGAGAGCCTGTTGATCTCCTGTCCACTCTGCCCTATGTTGGACCCTCAGTCA
 CACCAATGACCATCGTGGAAAACCTCAAGGAGTCTTCATTAAGGGCTTGTCAGTTCACCTCACAAGTCA
 GGAGGAGGATGCATTAGCCTCCTTTTTGAGGGTGAGACCAACAGGATTATAGCCTCCACACTATGAAC
 AAAAATCTTCCAGATCACACTGCATTTTACCATCTACTTAGAGGCCCATTCGCCGACCTATCAGAGG
 AAAAGTACATCACTTCCAAAATTAACCTGGTGGATCTGGCAGGCTCAGAGAGGCTGGGGAAGTCTGGGTC
 TGAGGGCCAAGTCTGAAGGAAGCCACCTACATCAACAAATCGCTCTCATTCTGGAGCAGGCCATCATT
 GCCCTTGGGGACCAAGCGGGACCATCCCCTTTCGGCAGTGAAGCTCACCCACGCTCTGAAGGACT
 CGTTAGGGGAACTGCAATATGGTCCTCGTGACAAACATCTATGGAGAAGCTGCCAGTTAGAAGAAAC
 GCTATCTTCACTGAGATTTGCCAGCAGGATGAAGCTAGTCACTGAGCCTGCCATCAATGAAAAGTAT
 GATGCTGAGAGAATGGTCAAGAACCTGGAGAAGGAAGTCACTCAAGCAGGAGCTGGCTATCCATG
 ACAGCTGACCAACCGCACCTTTGTGACCTATGACCCCATGGATGAAATCCAGATTGCTGAGATCAACTC
 CCAGGTGCGGAGGTACCTGGAGGGGACACTGGACGAGATCGACATAATCAGCCTTAGACAGATCAAGGAG
 GTGTTCAACCAAGTTCGGGTGGTCTGAGCCAACAGGAACAGGAAGTGGAGTCCACTTTGCGCAGGAAGT
 ACACCCTCATTGACAGGAATGACTTTGACGCCATTTCTGCTATCCAGAAGCGGGGCTTGTGGATGTTGA
 TGGCCACCTAGTGGGTGAGCCTGAAGGACAAAATTTGGACTCGGAGTCGCCCTTTCTCTACCAAACT
 GGAAGAAAGCCAAGTCCAAGAAGACATTCAAAGAGCCACTCAGCTCCTTGGCAAGAAAGGAAGGTGCCA
 GCAGCCCTGTGAATGGGAAGGACTTGGATTACGTTTCCACCTCCAAGACCCAGCTGGTCCCATCTCCAA
 AGATGGGGATGTCAAAGACATGCTTTGCGGGACCGGAACTTCCAGCATTGAGCCCCTTCCCTCAGAC
 TCCCCGAAGGAGGAATTACGCCAATTAGGCCGACACCCACCTCCAAACCAGTGGCCTTTGAGGAGT
 TTAAGAATGAGCAAGGTAGTGAGATCAACCGAATTTTCAAAGAAAACAAATCCATCTTGAATGAACGGAG
 GAAAAGGGCCAGCGAGACCACACAGCACATCAATGCCATCAAGCGGGAGATTGATGTGACCAAGGAGGCC
 CTGAATTTCCAGAAGTCACTACGGGAGAAGCAAGGCAAGTACGAAAACAAGGGGCTGATGATCATCGATG
 AGGAAGAATTCTGCTGATCCTCAAGCTCAAAGACCTCAAAGAAGCAGTACCGCAGCGAGTACCAGGACCT
 GCGTGACCTCAGGGCTGAGATCCAGTATTGCCAGCACCTAGTGGATCAGTGTGCCACCGCCTGCTCATG
 GAATTTGACATCTGGTACAATGAGTCCTTTGTCATCCCTGAGGACATGCAGATGGCACTGAAGCCAGGCG
 GCAGCATCCGGCCAGGCATGGTCCCTGTGAACAGGATTGTGTCTCTGGGAGAAGATGACCAGGACAAATT
 CAGCCAGCTGCAGCAGAGGGTGCTTCTGAGGGCCCTGATTCCATCTCCTTCTACAATGCCAAAGTCAAG
 ATAGAGCAGAAGCATAATTACTTGAAAACCATGATGGGCCTCCAGCAGGCACATAGAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226170 representing NM_001134878
 Red=Cloning site Green=Tags(s)


MGTRKKVHAFVRVKPTDDFAHEMIRYGDDKRSIDIHKKDIRRGVVNNQQTOWSFKLDGVLHDASQDLVY
 ETVAKDVVSQALDGYNGTIMCYGQTGAGKTYTMMGATENYKHRGILPRALQQVFRMIEERPTHAITVRVS
 YLEIYNESLFDLLSTLPYVGPSVTPMTIVENPQGVFIKGLSVHLTSQEEDAFSLLFEGETNRIIASHTMN
 KNSSRSHCIFTIYLEAHSRTLSEEKYITSKINLVDLAGSERLGKSGSEGQVLKEATYINKSLSFLEQAI
 ALGDQKRDPFRQCKLTHALKDSLGGNCNMVLVTNIYGEAAQLEETLSSLRFAARMKLVTEPAINEKY
 DAERMVKNLEKELALLKQELAIHDSL TNRTFVTYDPMDEIQIAEINSQVRRYLEGTLDEIDIIISLRQIKE
 VFNQFRVLSQQEQEVSTLRRKYTLIDRNDFAAISAIQKAGLVDVDGHLVGEPEGQNFGLGVAPFSTKP
 GKKAKSKKTFKEPLSSLARKEGASSPVNGKDLDYVSTSKTQLVPSSKDGVDKMLSRDRETSSIEPLPSD
 SPKEELRPIRPDTPPSKPVAFEEFKNEQGSEINRIFKENKSILNERRKRASETQHINAIKREIDVTKEA
 LNFQKSLREKQGYENKGLMIIDEEFLILKLKDLKKQYRSEYQDLRDLRAEIQYCQHLVDQCRHRLLM
 EFDIWNESFVIPEDMQMALKPGGSIRPGMVPVNRIVSLGEDDQDKFSQLQQRVLEPGPDSISFYNAKVK
 IEQKHNYLKTMMGLQQAHRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6318_g11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN:	NM_001134878
ORF Size:	2370 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001134878.2](#)

RefSeq ORF: 2373 bp

Locus ID: 64147


UniProt ID: [Q9HAQ2](#)

Cytogenetics: 3p21.31


Protein Families: Druggable Genome

MW: 89.8 kDa


Product images:



Circular map for RC226170



Western blot validation of overexpression lysate (Cat# [LY427510]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226170 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified KIF9 protein (Cat# [TP326170]). The protein was produced from HEK293T cells transfected with KIF9 cDNA clone (Cat# RC226170) using MegaTran 2.0 (Cat# [TT210002]).