

Product datasheet for RC226202L4

CCHCR1 (NM_001105563) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	mGFP
Symbol:	CCHCR1
Synonyms:	C6orf18; HCR; pg8; SBP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226202).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.

ACCN:	NM_001105563
ORF Size:	2505 bp



OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001105563.1 , NP_001099033.1
RefSeq ORF:	2508 bp
Locus ID:	54535
UniProt ID:	Q8TD31
Cytogenetics:	6p21.33
MW:	95 kDa
Gene Summary:	This gene encodes a protein with five coiled-coil alpha-helical rod domains that is thought to act as a regulator of mRNA metabolism through its interaction with mRNA-decapping protein 4. It localizes to P-bodies, the site of mRNA metabolism, with an N-terminus that is required for this subcellular localization, suggesting it is a P-body component. Naturally occurring mutations in this gene are associated with psoriasis. [provided by RefSeq, May 2017]

Product images:

Circular map for RC226202L4