

Product datasheet for **RC226219L3V**

FAM123C (AMER3) (NM_001105193) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM123C (AMER3) (NM_001105193) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM123C
Synonyms:	FAM123C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001105193
ORF Size:	2583 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226219).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001105193.1 , NP_001098663.1
RefSeq ORF:	2586 bp
Locus ID:	205147
UniProt ID:	Q8N944
Cytogenetics:	2q21.1
MW:	90.3 kDa
Gene Summary:	Regulator of the canonical Wnt signaling pathway. Acts by specifically binding phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5)P2), translocating to the cell membrane (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)