

Product datasheet for **RC226225L1V**

CCHCR1 (NM_001105564) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	CCHCR1
Synonyms:	C6orf18; HCR; pg8; SBP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001105564
ORF Size:	2613 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226225).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001105564.1 , NP_001099034.1
RefSeq Size:	2974 bp
RefSeq ORF:	2616 bp
Locus ID:	54535
UniProt ID:	Q8TD31
Cytogenetics:	6p21.33



MW: 99.5 kDa

Gene Summary: This gene encodes a protein with five coiled-coil alpha-helical rod domains that is thought to act as a regulator of mRNA metabolism through its interaction with mRNA-decapping protein 4. It localizes to P-bodies, the site of mRNA metabolism, with an N-terminus that is required for this subcellular localization, suggesting it is a P-body component. Naturally occurring mutations in this gene are associated with psoriasis. [provided by RefSeq, May 2017]