

Product datasheet for **RC226231L4V**

KIAA1310 (KANSL3) (NM_001115016) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KIAA1310 (KANSL3) (NM_001115016) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIAA1310
Synonyms:	KIAA1310; NSL3; Rcd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001115016
ORF Size:	2634 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226231).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001115016.2
RefSeq Size:	5240 bp
RefSeq ORF:	2637 bp
Locus ID:	55683
UniProt ID:	Q9P2N6
Cytogenetics:	2q11.2
MW:	93.2 kDa



[View online »](#)

Gene Summary:

As part of the NSL complex it is involved in acetylation of nucleosomal histone H4 on several lysine residues and therefore may be involved in the regulation of transcription.
[UniProtKB/Swiss-Prot Function]