

Product datasheet for **RC226241L2V**

CTNNA3 (NM_001127384) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CTNNA3 (NM_001127384) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CTNNA3 |
| Synonyms: | ARVD13; VR22 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001127384 |
| ORF Size: | 2685 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226241). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001127384.1 , NP_001120856.1 |
| RefSeq ORF: | 2688 bp |
| Locus ID: | 29119 |
| UniProt ID: | Q9UI47 |
| Cytogenetics: | 10q21.3 |
| Protein Pathways: | Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Endometrial cancer, Leukocyte transendothelial migration, Pathways in cancer, Tight junction |
| MW: | 99.6 kDa |



[View online »](#)

Gene Summary:

This gene encodes a protein that belongs to the vinculin/alpha-catenin family. The encoded protein plays a role in cell-cell adhesion in muscle cells. Mutations in this gene are associated with arrhythmogenic right ventricular dysplasia, familial 13. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014]