

Product datasheet for **RC226242L4V**

CASK (NM_001126055) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CASK (NM_001126055) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CASK
Synonyms:	CAGH39; CAMGUK; CMG; FGS4; hCASK; LIN2; MICPCH; MRXSNA; TNRC8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001126055
ORF Size:	2691 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226242).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001126055.2 , NP_001119527.1
RefSeq ORF:	2694 bp
Locus ID:	8573
UniProt ID:	O14936
Cytogenetics:	Xp11.4
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Tight junction
MW:	101.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a calcium/calmodulin-dependent serine protein kinase. The encoded protein is a MAGUK (membrane-associated guanylate kinase) protein family member. These proteins are scaffold proteins and the encoded protein is located at synapses in the brain. Mutations in this gene are associated with FG syndrome 4, intellectual disability and microcephaly with pontine and cerebellar hypoplasia, and a form of X-linked intellectual disability. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]