

## Product datasheet for **RC226243**

### CDKN1A interacting zinc finger protein 1 (CIZ1) (NM\_001131016) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CDKN1A interacting zinc finger protein 1 (CIZ1) (NM_001131016) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CDKN1A interacting zinc finger protein 1
Synonyms:	LSFR1; NP94; ZNF356
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RC226243 representing NM\_001131016  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC


ATGTTTCAGCCAGCAGCAGCAGCAGCAGCTCCAGCAACAGCAGCAGCAGCTCCAGCAGTTACAGCAGCAGC  
AGCTCCAGCAGCAGCAATTGCAGCAGCAGCAGTTACTGCAGCTCCAGCAGCTGCTCCAGCAGTCCCACC  
ACAGGCCCCGTTGCCATGGCTGTCAGCCGGGGCTCCCCCGCAGCAGCCACAGCAGCCGTTCTGAAT  
CTCCAGGGCACCAACTCAGCCTCCCTCCTCAACGGCTCCATGCTGCAGAGAGCTTTGCTTTTACAGCAGT  
TGCAAGGACTGGACCAGTTTGAATGCCACCAGCCACGTATGACACTGCCGGTCTCACCATGCCACAGC  
AACACTGGGTAACCTCCGAGGCTATGGCATGGCATCCCCAGGCCCTCGCAGCCCCCAGCCTCACACCCCCA  
CAACTGGCCACTCCAAATTTGCAACAGTTCTTTCCCCAGGCCACTCGCCAGTCTTGTGTTGGACCTCCTC  
CTGTTGGGGTCCCATGAACCTTCCCAGTTCAACCTTTTCAGGACGGAACCCCCAGAAACAGGCCGGAC  
CTCCTCCTTACCACCCCAATCGAAAGGATTCTTCTTCTCAGACAATGCCTGTGGAAGACAAGTCAGAC  
CCCCAGAGGGGTCTGAGGAAGCCGAGAGCCCCGGATGGACACACCAGAAGACCAAGATTTACCGCCT  
GCCCAGAGGACATCGCCAAGGAAAAACGCACTCCAGCACCTGAGCCTGAGCCTGTGAGGCGTCCGAGCT  
GCCAGCAAAGAGATTGAGGAGCTCAGAAGAGCCCACAGAGAAGGAACCTCCAGGGCAGTTACAGGTGAAG  
GCCCAGCCGAGGCCCGGATGACAGTACCGAAACAGACACAGACACCAGACCTGCTGCCTGAGGCCCTGG  
AAGCCCAAGTGTGCCACGATTCAGCCACGGGTCTGCAGGTCCAGGCCAGGTGCAGTCACAGACTCA  
GCCCGGATACCATCCACAGACACCCAGGTGCAGCCAAAGCTGCAGAAGCAGGCCAAACACAGACCTCT  
CCAGAGCACTTAGTGTGCAACAGAAGCAGGTGCAGCCACAGCTGCAGCAGGAGGCAGCCACAGAAGC  
AGGTGCAGCCACAGGTACAGCCACAGGCACATTCACAGGGCCCAAGGCAGGTGCAGTCAGCAGGAGGC  
AGAGCCGCTGAAGCAGGTGCAGCCACAGGTGCAGCCCCAGGCACATTCACAGCCCCAAGGCAGGTGCAG  
CTGCAGCTGCAGAAGCAGGTCCAGACACAGACATATCCACAGGTCCACACACAGGCACAGCCAAGCGTCC  
AGCCACAGGAGCATCCTCCAGCGCAGGTGTGAGTACAGCCACCAGAGCAGACCCATGAGCAGCCTCACAC  
CCAGCCGAGGTGTGTTGCTGGCTCCAGAGCAAACACCAGTTGTGGTTCATGTCTGCGGGCTGGAGATG  
CCACCTGATGCAGTAGAAGCTGGTGGAGGCATGGAAAAGACCTTGCCAGAGCCTGTGGGCACCCAAGTCA  
GCATGGAAGAGATTGAGAATGAGTCGGCTGTGGCCTAGATGTGGGAGAATGTGAAAACAGAGCGAGAGA  
GATGCCAGGGGTATGGGGCGCCGGGGCTCCCTGAAGGTACCATTCTGCAGAGCAGTGCAGCCGGGCC  
TTTAGCACTGTACCCCTGACACCTGTCCCCGCCCAAGTACTCGTCTCCTCCACCCCTGCGGCTACCA  
GCACTCCCTTAAGCAGGCCCTCCAGTTCTTCTGCTACATCTGCAAGGCCAGTGTCCAGCCAGCAGGA  
GTTCCAGGACCACATGTGCGGAGCCTCAGCACAGCAGCGCTAGGGGAGATCCAGCACATGAGCCAAGCC  
TGCCCTCCTGTCCCTGCTGCCCGTCCCCGGGACGTCTGGAGACAGAGGATGAGGAGCCTCCACCAAGGC  
GCTGGTGAACACCTGCCAGCTCTACTACATGGGGACCTGATCCAACACCCGAGGACACAGGACCACAA  
GATTGCCAAACAATCCTTGCGACCCTTCTGCACCGTTTGAACCGCTACTTCAAACCCCTCGCAAGTTT  
GTGGAGCACGTGAAGTCCCAGGGGCATAAAGGACAAAGCCAAGGAGCTGAAGTCGCTTGAGAAAGAAATTG  
CTGGCCAAAGATGAGGACCACTTCATTACAGTGGACGCTGTGGTTGCTTCGAGGGTGATGAAGAAGAGGA  
AGAGGATGATGAGGATGAAGAAGAGATCGAGGTTGAGGAGGAACCTGCAAGCAGGTGAGGTCCAGAGAT  
ATATCCAGAGAGGAGTGAAGGGCTCGGAGACCTACAGCCCAATACTGCATATGGTGTGACTTCTCTGG  
TGCCCGTATGGGCTATATCTGCCGCATCTGCCACAAGTTCTATCACAGCAACTCAGGGGCACAGCTCTC  
CCACTGCAAGTCCCTGGGCCACTTTGAGAACCTGCAGAAAACAAGGCGGCCAAGAACCACCCAGCCCCACC  
ACCCGACCTGTGAGCCGCGGTGCGCAATCAACGCCCGGAACGCTTTGACAGCCCTGTTACCTCCAGCG  
GCCGCCACCCTCCAGCCCAACACCCAGGACAAAACACCCAGCAAGGTGACGGCTCGACCCTCCAGCC  
CCCACTACCTCGGCGCTCAACCCGCTCAAACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA




<b>ORF Size:</b>	2694 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001131016.2</a>
<b>RefSeq ORF:</b>	2697 bp
<b>Locus ID:</b>	25792
<b>UniProt ID:</b>	<a href="#">Q9ULV3</a>
<b>Cytogenetics:</b>	9q34.11
<b>MW:</b>	99.9 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a zinc finger DNA binding protein that interacts with CIP1, part of a complex with cyclin E. The encoded protein may regulate the cellular localization of CIP1. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]

Product images:



Circular map for RC226243



Western blot validation of overexpression lysate (Cat# [LY427357]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226243 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).