

Product datasheet for **RC226244L4V**

CASK (NM_001126054) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CASK (NM_001126054) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CASK |
| Synonyms: | CAGH39; CAMGUK; CMG; FGS4; hCASK; LIN2; MICPCH; MRXSNA; TNRC8 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001126054 |
| ORF Size: | 2694 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226244). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001126054.2 , NP_001119526.1 |
| RefSeq ORF: | 2697 bp |
| Locus ID: | 8573 |
| UniProt ID: | O14936 |
| Cytogenetics: | Xp11.4 |
| Protein Families: | Druggable Genome, Protein Kinase |
| Protein Pathways: | Tight junction |
| MW: | 102 kDa |



[View online »](#)

Gene Summary:

This gene encodes a calcium/calmodulin-dependent serine protein kinase. The encoded protein is a MAGUK (membrane-associated guanylate kinase) protein family member. These proteins are scaffold proteins and the encoded protein is located at synapses in the brain. Mutations in this gene are associated with FG syndrome 4, intellectual disability and microcephaly with pontine and cerebellar hypoplasia, and a form of X-linked intellectual disability. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]