

## Product datasheet for **RC226257L2V**

### SLFN14 (NM\_001129820) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLFN14 (NM_001129820) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLFN14
Synonyms:	BDPLT20
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001129820
ORF Size:	2736 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226257).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001129820.1</a>
RefSeq Size:	2882 bp
RefSeq ORF:	2739 bp
Locus ID:	342618
UniProt ID:	<a href="#">P0C7P3</a>
Cytogenetics:	17q12
MW:	103.9 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene plays an important role in platelet formation and function. Defects in this gene are a cause of thrombocytopenia with excessive bleeding. [provided by RefSeq, Jul 2016]