

Product datasheet for **RC226263L3V**

TRPC3 (NM_001130698) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRPC3 (NM_001130698) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRPC3
Synonyms:	SCA41; TRP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130698
ORF Size:	2763 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226263).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130698.1 , NP_001124170.1
RefSeq ORF:	2766 bp
Locus ID:	7222
UniProt ID:	Q13507
Cytogenetics:	4q27
Protein Families:	Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane
MW:	105.4 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a membrane protein that can form a non-selective channel permeable to calcium and other cations. The encoded protein appears to be induced to form channels by a receptor tyrosine kinase-activated phosphatidylinositol second messenger system and also by depletion of intracellular calcium stores. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]