

## Product datasheet for **RC226265L4V**

### NCOA7 (NM\_001122842) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NCOA7 (NM_001122842) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NCOA7
Synonyms:	dj187J11.3; ERAP140; ESNA1; Nbla00052; Nbla10993; NCOA7-AS; TLDC4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001122842
ORF Size:	2793 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226265).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001122842.2</a>
RefSeq Size:	6424 bp
RefSeq ORF:	2796 bp
Locus ID:	135112
UniProt ID:	<a href="#">Q8NI08</a>
Cytogenetics:	6q22.31-q22.32
Protein Families:	Druggable Genome
MW:	104.8 kDa


[View online »](#)

**Gene Summary:**

Enhances the transcriptional activities of several nuclear receptors. Involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA.  
[UniProtKB/Swiss-Prot Function]