



Product datasheet for RC226268L4

FPGT-TNNI3K (NM_001112808) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FPGT-TNNI3K (NM_001112808) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	FPGT-TNNI3K
Synonyms:	CARK; TNNI3K
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226268).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.

ACCN:	NM_001112808
ORF Size:	2847 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001112808.2, NP_001106279.2</u>
RefSeq Size:	3328 bp
RefSeq ORF:	2811 bp
Locus ID:	100526835
Cytogenetics:	1p31.1
MW:	105.6 kDa
Gene Summary:	This locus represents naturally occurring read-through transcription from the neighboring fucose-1-phosphate guanylyltransferase (FPGT) and TNNI3 interacting kinase (TNNI3K) genes. Alternative splicing results in multiple transcript variants that are composed of in-frame exons from each individual gene. [provided by RefSeq, Dec 2010]



This product is to be used for laboratory only. Not for diagnostic or therapeutic use.
©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US