

Product datasheet for **RC226285L4V**

SYNPO2L (NM_001114133) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SYNPO2L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001114133
ORF Size:	2931 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226285).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001114133.1
RefSeq ORF:	2934 bp
Locus ID:	79933
UniProt ID:	Q9H987
Cytogenetics:	10q22.2
MW:	102.3 kDa
Gene Summary:	Actin-associated protein that may play a role in modulating actin-based shape. [UniProtKB/Swiss-Prot Function]



