

Product datasheet for **RC226289L3V**

PIGG (NM_001127178) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	PIGG
Synonyms:	GPI7; LAS21; MRT53; PRO4405; RLGS1930
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001127178
ORF Size:	2949 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226289).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001127178.1
RefSeq Size:	3266 bp
RefSeq ORF:	2952 bp
Locus ID:	54872
UniProt ID:	Q5H8A4
Cytogenetics:	4p16.3



Protein Families: Transmembrane

Protein Pathways: Glycosylphosphatidylinositol(GPI)-anchor biosynthesis

MW: 108.6 kDa

Gene Summary: This gene encodes an enzyme involved in glycosylphosphatidylinositol-anchor biosynthesis. The encoded protein, which is localized to the endoplasmic reticulum, is involved in transferring ethanolamine phosphate to mannose 2 of glycosylphosphatidylinositol species H7 to form species H8. Allelic variants of this gene have been associated with intellectual disability, hypotonia, and early-onset seizures. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]