

## Product datasheet for **RC226299L4V**

### IPO11 (NM\_001134779) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IPO11 (NM_001134779) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IPO11
Synonyms:	RanBP11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001134779
ORF Size:	3045 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226299).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001134779.1</a> , <a href="#">NP_001128251.1</a>
RefSeq ORF:	3048 bp
Locus ID:	51194
UniProt ID:	<a href="#">Q9UI26</a>
Cytogenetics:	5q12.1
MW:	116.8 kDa



[View online »](#)

**Gene Summary:**

Importins, including IPO11, are a members of the karyopherin/importin-beta family of transport receptors (see KPNB1; 602738) that mediate nucleocytoplasmic transport of protein and RNA cargoes (Plafker and Macara, 2000 [PubMed 11032817]).[supplied by OMIM, Sep 2008]