

Product datasheet for **RC226304L4V**

SLC12A3 (NM_001126107) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC12A3 (NM_001126107) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC12A3
Synonyms:	NCC; NCCT; TSC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001126107
ORF Size:	3087 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226304).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001126107.1 , NP_001119579.1
RefSeq ORF:	3090 bp
Locus ID:	6559
UniProt ID:	P55017
Cytogenetics:	16q13
Protein Families:	Druggable Genome, Transmembrane
MW:	113.6 kDa


[View online »](#)

Gene Summary:

This gene encodes a renal thiazide-sensitive sodium-chloride cotransporter that is important for electrolyte homeostasis. This cotransporter mediates sodium and chloride reabsorption in the distal convoluted tubule. Mutations in this gene cause Gitelman syndrome, a disease similar to Bartter's syndrome, that is characterized by hypokalemic alkalosis combined with hypomagnesemia, low urinary calcium, and increased renin activity associated with normal blood pressure. This cotransporter is the target for thiazide diuretics that are used for treating high blood pressure. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]