

Product datasheet for **RC226307L3V**

FAM160A1 (NM_001109977) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	FAM160A1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001109977
ORF Size:	3120 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001109977.1
RefSeq ORF:	3123 bp
Locus ID:	729830
UniProt ID:	Q05DH4
Cytogenetics:	4q31.3
MW:	116.4 kDa



Gene Summary:

Probable component of the FTS/Hook/FHIP complex (FHF complex) (PubMed:32073997). FHF complex promotes the distribution of AP-4 complex to the perinuclear area of the cell (PubMed:32073997).[UniProtKB/Swiss-Prot Function]