

## Product datasheet for **RC226313L3V**

### LIMCH1 (NM\_001112717) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | LIMCH1 (NM_001112717) Human Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                   |
| Symbol:                   | LIMCH1                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Synonyms:                 | LIMCH1A; LMO7B                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                                                                                                                                                                                                                                                                                                                                                                                                               |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_001112717                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| ORF Size:                 | 3171 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC226313).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_001112717.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                     |
| RefSeq ORF:               | 3174 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:                 | 22998                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| UniProt ID:               | <a href="#">Q9UPQ0</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:             | 4p13                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| MW:                       | 118.9 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                          |


[View online »](#)

**Gene Summary:**

Actin stress fibers-associated protein that activates non-muscle myosin IIa. Activates the non-muscle myosin IIa complex by promoting the phosphorylation of its regulatory subunit MRLC/MYL9. Through the activation of non-muscle myosin IIa, positively regulates actin stress fibers assembly and stabilizes focal adhesions. It therefore negatively regulates cell spreading and cell migration.[UniProtKB/Swiss-Prot Function]