

## Product datasheet for **RC226360L3V**

### **FNDC3B (NM\_001135095) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FNDC3B (NM_001135095) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FNDC3B
Synonyms:	FAD104; PRO4979; YVTM2421
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135095
ORF Size:	3612 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226360).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001135095.1</a> , <a href="#">NP_001128567.1</a>
RefSeq ORF:	3615 bp
Locus ID:	64778
UniProt ID:	<a href="#">Q53EP0</a>
Cytogenetics:	3q26.31
Protein Families:	Druggable Genome, Transmembrane
MW:	132.7 kDa
Gene Summary:	May be a positive regulator of adipogenesis.[UniProtKB/Swiss-Prot Function]


[View online »](#)