

## Product datasheet for **RC226364L4V**

### **RBM20 (NM\_001134363) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RBM20 (NM_001134363) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBM20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001134363
ORF Size:	3681 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226364).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001134363.1</a> , <a href="#">NP_001127835.1</a>
RefSeq ORF:	3684 bp
Locus ID:	282996
UniProt ID:	<a href="#">Q5T481</a>
Cytogenetics:	10q25.2
MW:	134.2 kDa
Gene Summary:	This gene encodes a protein that binds RNA and regulates splicing. Mutations in this gene have been associated with familial dilated cardiomyopathy. [provided by RefSeq, Apr 2014]



[View online »](#)