

Product datasheet for RC226373L4

DOCK9 (NM_001130050) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DOCK9 (NM_001130050) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	DOCK9
Synonyms:	ZIZ1; ZIZIMIN1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226373).
Restriction Sites:	SgfI-RsrII
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001130050
ORF Size:	3759 bp



[View online »](#)

