

Product datasheet for **RC226390L2V**

ROBO2 (NM_001128929) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | ROBO2 (NM_001128929) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ROBO2 |
| Synonyms: | SAX3 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001128929 |
| ORF Size: | 4182 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226390). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001128929.1 , NP_001122401.1 |
| RefSeq ORF: | 4185 bp |
| Locus ID: | 6092 |
| UniProt ID: | Q9HCK4 |
| Cytogenetics: | 3p12.3 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Axon guidance |
| MW: | 152.9 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the ROBO family, part of the immunoglobulin superfamily of proteins that are highly conserved from fly to human. The encoded protein is a transmembrane receptor for the slit homolog 2 protein and functions in axon guidance and cell migration. Mutations in this gene are associated with vesicoureteral reflux, characterized by the backward flow of urine from the bladder into the ureters or the kidney. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014]